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- 1. 河原玲子：日本臨牀、47(増刊)1033-1035(1989)
- 2. Koenig RJ, Peterson CM, Kilo C, Cermami A, Williamson JR. : Hemoglobin A1c as an indicator of the degree of glucose intolerance in diabetes. *Diabetes*, 25(3):230-2(1976)
- 3. American Diabetes Association. Implications of the United Kingdom Prospective Diabetes Study(Position Statement). *Diabetes Care*, 22:27-31(1999)
- 4. American Diabetes Association. Tests of Glycemia in diabetes (Position Statement). *Diabetes Care*, 21:77-79(1999)
- 5. Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*, 329:977-86(1993)
- 6. Diabetes Control and Complications Trial Research Group. The Absence of a glycemic threshold for the development of long-term complications: The perspective of the diabetes control and complications trial. *Diabetes*, 45:1289-98(1996)
- 7. Krishnamurti U, Steffes MW. Glycohemoglobin: A primary predictor of the development or reversal of complications of diabetes mellitus. *Clin Chem*, 47:1157-1165(2001)
- 8. Kashiwagi A, et al. International clinical harmonization of glycated hemoglobin in Japan: From Japan Diabetes Society to national glycohemoglobin standardization program values. *Diabetology International*. 3:8-10.
- 9. Kashiwagi A, et al. International clinical harmonization of glycated hemoglobin in Japan: From Japan Diabetes Society to national glycohemoglobin standardization program values. *Journal of Diabetes Investigation*. 3:39-40.
- 10. 島 健二ほか：グリコヘモグロビンの標準化に関する委員会の中間報告、糖尿病、37：233-241(1994)

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