HbA1c

081248/03-2007

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1. A. Finke, U. Kobold, W. Hoelzel, C. Weykamp, K. Miedema, J-O. Jeppsson. Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood. Clin Chem Lab Med; 40(1): 78-89, 2002.
2. Goldstein DE, et al, Is glycohemoglobin testing useful in diabetes mellitus: Lesson from the diabetes control and complication trial. Clin Chem 1994; 40: 1637-1640.
3. The Diabetes Control and Complications Trial Research Group: The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med 329: 977-986; 1993.
4. Intensive blood glucose control with sulphonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes(UKPDS33). UK Prospective Diabetes Study(UKPDS) Group. Lancet 325: 837-853; 1998.
5. Burtis C, et al, Text book of Clinical Chemistry, 3rd ed. Philadelphia, PA; WB Saunders Co; 791, 1999.
6. Calam RR. Specimen Processing Separator Gels: An Update. J Clin Immunoassay. 11: 86-90; 1988.
7. Kashiwagi A., et al. International clinical harmonization of glycated hemoglobin in Japan: From Japan Diabetes Society to national glycohemoglobin standardization program values. Diabetology International. 3:8-10.
8. Kashiwagi A., et al. International clinical harmonization of glycated hemoglobin in Japan: From Japan Diabetes Society to national glycohemoglobin standardization program values. Journal of Diabetes Investigation. 3:39-40.
9. <http://diabetes.niddk.nih.gov/dm/pubs/hemovari-A1c/index.aspx#hemoglobinopathies>. Accessed on April 28, 2014.
10. <http://www.ngsp.org/factors.asp>. Accessed on April 28, 2014.
11. Bry L, Chen PC, Sacks DB. Effects of hemoglobin variants and chemically modified derivatives on assays for glycated hemoglobin. Clin Chem. 47: 153-163; 2001.
12. International Association of Diabetes and Pregnancy Study Groups Recommendations on the Diagnosis and Classification of Hyperglycemia in Pregnancy. Diabetes Care. 33:No 3; 2010.
13. International Expert Committee Report on the Role of the A1c Assay in the Diagnosis of Diabetes. Diabetes Care. 32(7):1327-1334; 2009.
14. Standards of Medical Care in Diabetes. Diabetes Care. 38:Supplement 1; S8-S16, 2015.

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