PSA II Reag

059339/04-2012

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1. 平成24年 人口動態統計 (厚生労働省大臣官房統計情報部)
2. Wang MC, Valenzuela LA, Murphy GP and Chu TM. Purification of a human prostate specific antigen. Invest Urol. 17:159-163; 1979.
3. Wang MC, Papsidero LD, Kuriyama M et al. Prostate antigen: a new potential marker for prostatic cancer. The Prostate. 2:89-96; 1981.
4. Papsidero LD, Wang MC, Valenzuela LA et al. A Prostate Antigen in Sera of Prostatic Cancer Patients. Cancer Res. 40:2428-2432; 1980.
5. Oesterling JE. Prostate specific antigen: a critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. J Urol. 145:907-923; 1991.
6. Stamey TA, Kabalin JN, McNeal JE et al. Prostate specific antigen in the diagnosis and treatment of adenocarcinoma of the prostate. II. Radical prostatectomy treated patients. J Urol. 141:1076-1083; 1989.
7. Chybowski FM, Larson Keller JJ, Bergstrahl EJ and Oestriling JE. Predicting radionuclide bone scan findings in patients with newly diagnosed, untreated prostate cancer: prostate specific antigen is superior to all other clinical parameters. J Urol. 145:313-318; 1991.
8. Stamey TA, Kabalin JN, Ferrari M. Prostate specific antigen in the diagnosis and treatment of adenocarcinoma of the prostate. III. Radiation treated patients. J Urol. 141:1084-1087; 1989.
9. Stamey TA, Kabalin JN, Ferrari M, Yang N. Prostate specific antigen in the diagnosis and treatment of adenocarcinoma of the prostate. IV. Anti-androgen treated patients. J Urol. 141:1088-1090; 1989.
10. Summers M et al. Luminogenic Reagent Using 3-Chloro 4-Hydroxy Acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. Clin Chem. 41:S73; 1995.
11. Calam RR. Specimen Processing Separator Gels: An Update. J Clin Immunoassay. 11:86-90; 1988
12. Levinson, SS.: The Nature of Heterophilic Antibodies and Their Role in Immunoassay Interference. J Clin Immunoassay 15: 108-115 (1992).
13. Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry; Approved Guideline, Second Edition. CLSI document EP07-A2. Wayne, PA: CLSI; 2005.
14. Livaniou E, et.al.: Serum biotin levels in patients undergoing chronic hemodialysis. Nephron. 1987; 46(3):331-2.
15. Paul Grimsey et.al.: Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. International Journal of Pharmacokinetics. 2017; 2 (4): 247-56.
16. European Commission Scientific Committee on Food: Opinion of the Scientific Committee on Food on the Tolerable Upper Intake Level of Biotin (expressed on 26 September 2001). Rue de la Loi 200, B-1049 Bruxelles/Wetstraat 200, B-1049 Brussels-Belgium.
17. Mock and Mock: Serum concentrations of bisnorbiotin and biotin sulfoxide increase during both acute and chronic biotin supplementation J Lab Clin Med 1997;129:384-8.
18. Marcelo C. Batista: Biotin interference in immunoassays mimicking subclinical Graves' disease and hyperestrogenism: a case series. Clin Chem Lab Med. 2017 May 1;55(6):e99-e103.
19. Young DS. Effects of Drugs on Clinical Laboratory Tests. ed. 4. Washington, D.C.: AACC Press; 1995.
20. Friedman RB, Young DS. Effects of Disease on Clinical Laboratory Tests. ed. 3. Washington, D.C.: AACC Press; 1997.
21. Tryding N, Tufvesson C, Sonntag O (eds). Drug Effects in Clinical Chemistry. ed. 7. Stockholm: The National Corporation of Swedish Pharmacies, Pharmasoft AB, Swedish Society for Clinical Chemistry; 1996.
22. Guess, H. A. et al. : J. Urol., 155(1) : 3, 1996
23. Andriole GL et al : J Urol 175(5) : 1657-1662, 2006
24. Marks LS et al : J Urol 176(3) : 868-874, 2006
25. Passing HW, A New Biometrical Procedure of Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Biochem. 21:709-720 1983.

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