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- Fairbanks VF, Ziesmer SC, O'Brien PC. Methods for measuring plasma hemoglobin in micromolar concentration compared. *Clin Chem*. 1992;38 (1):132-140.
- Knowles BJ, Haigh WB, Michaud GC. A monoclonal antibody-based immunoassay for hemoglobin A1c. *Diabetes*. 1986;35 (supplement 94A).
- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.

6. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline-Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
7. Data on file at Siemens Healthcare Diagnostics.
8. Wu AHB, ed. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Saunders; 2006; 480-483.
9. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
10. Bry L, Chen PC, Sacks DB. Effects of hemoglobin variants and chemically modified derivatives on assays for glycohemoglobin. *Clin Chem*. 2001;47 (2) :153-163.
11. Little RR, Rohlfing CL, Hanson SE, Roberts WL. Effects of hemoglobin C and S traits on Tosoh G8 and Siemens Advia HbA1c assays. *Clin Chim Acta*. 2010;411 (9-10) :779-780.
12. Hoelzel W, Weykamp C, Jeppsson JO, et al. IFCC reference system for measurement of hemoglobin A1c in human blood and the national standardization schemes in the United States, Japan, and Sweden: A method-comparison study. *Clin Chem*. 2004;50 (1) :166-174.
13. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
14. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. Washington, DC: AACC Press; 2000.
15. Rohlfing CL, Wiedmeyer HM, Little RR, England JD, Tennill A, Goldstein DE. Defining the relationship between plasma glucose and HbA (1c) : Analysis of glucose profiles and HbA (1c) in the Diabetes Control and Complications Trial. *Diabetes Care*. 2002;25 (2) : 275-178.
16. The Diabetes Control and Complications Trial Research Group. The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. *N Engl J Med*. 1993;329 (14) :977-86.
17. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
18. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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