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1. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
2. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
3. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
4. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
5. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*. 3rd ed. Washington, DC: AACC Press; 2007:529.
6. Schreiber W, Rausch S, Lammers M. Choice of specimen for immunonephelometric protein assays. *Clin Chim Acta*. 2005;355:S407.
7. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
8. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
9. Tietz NW. *Clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:354-355.
10. a. Heidelberger M, Kendall FE. *J Exper Med*. 1935;61:563-591.
b. *J Exper Med*. 1935;62:467-483.
c. *J Exper Med*. 1935;63:697-720.
11. Libby RL. A new and rapid quantitative technic for the determination of potency of types I and II antipneumococcal serum. *J Immunol*. 1938;34:269-279.
12. Gitlin D, Edelhoeh H. A study of the reaction between human serum albumin and its homologous equine antibody through the medium of light scattering. *J Immunol*. 1951;66 (1) :67-77.
13. Goldberg RJ, Campbell DH. The light-scattering properties of an antigen-antibody reaction. *J Immunol*. 1951;66:79-86.
14. Johnson P, Ottewill RH. A light scattering study of diphtheria toxin-antitoxin interaction. *Discuss Faraday Soc*. 1954;18:327-337.
15. Marrack JR, Grant RA. The interaction of antigens and antibodies in low concentrations of salt. *Br J Exp Pathol*. 1953;34 (3) :263-272.
16. Marrack JR, Richards CB. Light-scattering studies of the formation of aggregates in mixtures of antigen and antibody. *Immunology*. 1971;20 (6) :1019-1040.
17. Ritchie RF. A simple, direct, and sensitive technique for measurement of specific protein in dilute solution. *J Lab Clin Med*. 1967;70 (3) :512-517.
18. Helsing K. Influence of polymers on the antigen-antibody reaction in a continuous flow system. In: *Automated Immunoprecipitin Reaction. Colloquium on AIP*. Tarrytown, NY: Technicon Inst. Corp.; 1972:17.
19. Helsing K, Laurent TC. The influence of dextran on the precipitin reaction. *Acta Chem Scand*. 1964;18:1303.
20. Helsing K. The effect of dextran on the reaction between iodine-125 labelled human serum albumin and gamma-G-globulin from rabbit anti-albumin sera. *Acta Chem Scand*. 1966;20 (5) :1251-1262.
21. Helsing K. The effect of hyaluronate, chondroitin sulphate and chondroitin sulphate-protein complex on the precipitin reaction. *Biochem J*. 1969;112 (4) :475-481.

22. Hellsing K. Polysaccharide-enhanced precipitin reactions with antigens of various sizes. *Biochem J.* 1969;114 (1) :145-149.
23. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- * * 24. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.
25. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
26. Data on file at Siemens Healthcare Diagnostics.
27. Young DS. *Effects of Drugs on Clinical Laboratory Tests.* 3rd ed. Washington, DC: AACC Press; 1990.

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