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1. Kaplan LA, Pesce AJ, eds. *Clinical Chemistry: Theory, Analysis and Correlation*. St. Louis, MO: Mosby Company; 1989.
2. Price CP, Hammond PM, Scawen MD. Evaluation of an enzymic procedure for the measurement of acetaminophen. *Clin Chem*. 1983;29 (2) :358-361.
3. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
4. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
5. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
6. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
7. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*. 3rd ed. Washington, DC: AACC Press, 2007;46.
8. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
9. Data on file at Siemens Healthcare Diagnostics.
10. Wu AHB, ed. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Saunders; 2006:1238.
11. Merck Manual. <https://www.merckmanuals.com/professional/injuries-poisoning/poisoning/acetaminophen-poisoning>.
12. Ameer B, Greenblatt DJ. Acetaminophen. *Ann Intern Med*. 1977;87 (2) :202-209.
13. Barker JD, de Carle DJ, Anuras S. Chronic excessive acetaminophen use and liver damage. *Ann Intern Med*. 1977;87 (3) :299-301.

14. Prescott LF. The nephrotoxicity and hepatotoxicity of antipyretic analgesics. *Br J Clin Pharmacol*. 1979;7 (5) :453-462.
15. Black M. Acetaminophen hepatotoxicity. *Gastroenterology*. 1980;78 (2) :382-392.
16. Ambre J, Alexander M. Liver toxicity after acetaminophen ingestion. Inadequacy of the dose estimate as an index of risk. *JAMA*. 1977;238 (6) :500-501.
17. Vale JA, Meredith TJ, Goulding R. Treatment of acetaminophen poisoning. The use of oral methionine. *Arch Intern Med*. 1981;141 (3 Spec No) :394-396.
18. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- *19. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.
20. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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