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- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Wu AHB. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Saunders; 2006:142.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Brown SA, Hutchinson R, Morrisett J, Boerwinkle E, et al. Plasma Lipid, Lipoprotein Cholesterol, and APOprotein Distribution in Selected US Communities. *Arteriosclerosis and Thrombosis*. 13:1139-1158 (1993).
- Tietz NW. *Clinical Guide to Laboratory Tests*. 3rd ed. Philadelphia, PA: WB Saunders Company; 1995:68-69.
- Alaupovic P. Structure and function of plasma lipoproteins with particular regard to hyperlipoproteinemias and atherosclerosis. *Ann Biol Clin*. 1980;38 (2) :83-93.
- Naito HK. The clinical significance of apolipoprotein measurements. *J Clin Immunoassay*. 1986;9:11-20.
- Burtis CA, Ashwood ER, Bruns DE, eds. *Tietz Fundamentals of Clinical Chemistry*. 6th ed. St. Louis, MO: Saunders; 2008.
- Riesen WF, Mordasini R, Salzmann C, Theler A, Gurtner HP. Apoproteins and lipids as discriminators of severity of coronary heart disease. *Atherosclerosis*. 1980;37 (1) :157-162.
- Sandkamp M. Apolipoprotein-Diagnostik: Klinische Relevanz. *Diagnose & Labor*. 1990;40:37-43.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.
- Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- Data on file at Siemens Healthcare Diagnostics.

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