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1. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
2. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
3. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
4. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
5. Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*. 3rd ed. Washington, DC: AACC Press, 2007;217.
6. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
7. Data on file at Siemens Healthcare Diagnostics.
8. Parant F, Bossu H, Gagnieu MC, Lardet G, Moulsmas M. Cross-reactivity assessment of carbamazepine-10,11-epoxide, oxcarbazepine, and 10-hydroxycarbamazepine in two automated carbamazepine immunoassays: PETINIA and EMIT 2000. *Ther Drug Monit*. 2003;25 (1) :41-45.
9. Gilman A, Goodman L, Rall TW, Murad F. *Goodman and Gilman's The Pharmacological Basis of Therapeutics*. 7th ed. New York, NY: MacMillan Publishing Co; 1985:458-459.
10. Burtis CA, Ashwood ER, eds. *Tietz Textbook of Clinical Chemistry*. 2nd ed. Philadelphia: WB Saunders Co.; 1994;2212.
11. Baselt RC. *Disposition of Toxic Drugs and Chemicals in Man*. 5th ed. Foster City, CA: Chemical Toxicology Institute; 2000:114.
12. Kaplan LA, Pesce AJ, eds. *Clinical Chemistry: Theory, Analysis, Correlation*. 2nd ed. St. Louis, MO: C.V. Mosby Company; 1989:844.
13. Shen S, Elin RJ, Soldin SJ. Characterization of cross-reactivity by carbamazepine-10,11-epoxide with carbamazepine assays. *Clin Biochem*. 2001;34 (2) :157-158.
14. Bertilsson L, Tomson T. Clinical pharmacokinetics and pharmacological effects of carbamazepine and carbamazepine-10,11-epoxide: An update. *Clin Pharmacokinet*. 1986;11 (3) :177-198.
15. Hundt HK, Aucamp AK, Müller FO, Potgieter MA. Carbamazepine and its major metabolites in plasma: A summary of eight years of therapeutic drug monitoring. *Ther Drug Monit*. 1983;5 (4) :427-35.
16. Taylor WJ, Finn AL, eds. *Individualizing Drug Therapy: Practical Applications of Drug Monitoring*. Vol 2. New York: Gross, Townsend, Frank Inc; 1981:1-25.
17. Burtis CA, Ashwood ER, Bruns DE, eds. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics*. 4th ed. St. Louis, MO: Saunders; 2006:1249,2305.

18. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

\* \* 19. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.

20. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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