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- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Schreiber W, Rausch S, Lammers M. Choice of specimen for immunonephelometric protein assays. *Clin Chim Acta*. 2005;355:S407.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*. Washington, DC: AACC Press, 2007:537.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Wu AHB, ed. *Tietz Clinical Guide to Laboratory Tests*. 4th ed. St. Louis, MO: Saunders; 2006:608.
- a. Heidelberger M, Kendall FE. *J Exper Med*. 1935;61:563-591.  
b. *J Exper Med*. 1935;62:467-483.  
c. *J Exper Med*. 1935;63:697-720.
- Libby RL. New and rapid quantitation technique for the determination of potency of types I and II antipneumococcal serum. *J Immunol*. 1938;34:269.
- Gitlin D, Edelhock H. A study of the reaction between human serum and its homologous equine antibody through the medium of light scattering. *J Immunol*. 1951;66:67-77.
- Goldberg RJ, Campbell DH. The light scattering properties of an antigen-antibody reaction. *J Immunol*. 1951;66:70.
- Johnson P, Ottewill RH. A light scattering study of diphtheria toxin-antitoxin interaction. *Disc Faraday Soc*. 1954;18:327.
- Marrack JR, Grant RA. The interaction of antigen and antibody in low concentrations of salt. *Brit J Exp Path*. 1953;34:263.
- Marrack JR, Richards CB. Light scattering studies of the formation of aggregates in mixtures of antigen and antibody. *J Immunol* 1971;20:1019-1038.
- Ritchie RF. A simple, direct and sensitive technique for measurement of specific protein in dilute solution. *J Lab Clin Med*. 1967;70:512-517.
- Helsing K. Influence of polymers on the antigen-antibody reaction in a continuous flow system. In: *Automated Immunoprecipitin Reaction. Colloquium on AIP*. Tarrytown, NY: Technicon Inst. Corp.; 1972:17.
- Helsing K, Laurent TC. The influence of dextran on the precipitin reaction. *Acta Chem Scand*. 1964;18:1303.
- Helsing K. The effect of dextran on the reaction between iodine-125 labeled human serum albumin and gamma G-globulin from rabbit anti-albumin sera. *Acta Chem Scand*. 1966;20:1251.
- Helsing K. The effect of hyaluronate, chondroitin sulfate and chondroitin sulphate-protein complex on the precipitin reaction. *Biochem J*. 1969;112:475.
- Helsing K. Polysaccharide-enhanced precipitin reactions with antigens of various sizes. *Biochem J*. 1969;114:145.

23. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- \* \* 24. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.
25. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
26. Data on file at Siemens Healthcare Diagnostics.
27. Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. Washington, DC: AACC Press; 2000.

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