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- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. CLSI Document H18-A3.
- Data on file at Siemens Healthcare Diagnostics.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Burtis CA, Ashwood ER, Bruns DE. *Tietz Textbook of Clinical Chemistry and Molecular Biology*. 4th ed. St. Louis, MO: Saunders; 2006:1253,2315.
- Finn AL, Taylor WJ. *Individualizing Drug Therapy, Practical Applications of Drug Monitoring*. New York: Gross, Townsend, Frank, Inc.; 1981:87-108.
- Private Communication. Lilly Research Laboratories. A Division of Eli Lilly and Company, Indianapolis, Indiana, February 1993.
- Freeman CD, Quintiliani R, Nightingale CH. Vancomycin drug monitoring: Is it necessary? *Ann Pharmacother*. 1993;27 (5) :594-598.
- Fitzsimmons WE, Postelnick MJ, Tortorice PV. Survey of vancomycin monitoring guidelines in Illinois hospitals. *Drug Intel Clin Pharm*. 1988;22 (7-8) :598-600.
- Healy DP, Polk RE, Garson ML, Rock DT, Comstock TJ. Comparison of steady-state pharmacokinetics of two dosage regimens of vancomycin in normal volunteers. *Antimicrob Agents Chemother*. 1987;31 (3) :393-397.
- Lake KD, Peterson CD. Evaluation of a method for initiating vancomycin therapy: Experience in 205 patients. *Pharmacotherapy*. 1988;8 (5) :284-286.
- Blouin RA, Bauer LA, Miller DD, Record KE, Griffen WO Jr. Vancomycin pharmacokinetics in normal and morbidly obese subjects. *Antimicrob Agents Chemother*. 1982;21 (4) :575-580.
- Cheung RP, DiPiro JT. Vancomycin: An update. *Pharmacotherapy*. 1986;6 (4) :153-169.
- Polk RE, Espinel-Ingroff A, Lockridge R. In vitro evaluation of a vancomycin radioimmunoassay and observations on vancomycin pharmacokinetics in dialysis patients. *Drug Intel Clin Pharm*. 1981;15 (1) :15-20.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests*. 5th ed. Washington, DC: AACC Press; 2000.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.

21. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- \*22. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP07.

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