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- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Clinical and Laboratory Standards Institute. *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2000. CLSI Document C28-A2.
- Carter HB, Albertsen PC, Barry MJ; American Urological Association. Early detection of prostate cancer: AUA Guideline. *J Urol*. 2013;190 (2) :419-426.
- Mottet N, Bellmunt J, Briers E, et al; European Association of Urology. EAU Guidelines- *Prostate Cancer*. Arnhem, The Netherlands: European Association of Urology; 2017. <https://uroweb.org/guideline/prostate-cancer>.
- Catalona WJ, Smith DS, Ratliff TL, et al. Measurement of prostate-specific antigen in serum as a screening test for prostate cancer. *N Engl J Med*. 1991;324 (17) :1156-1161.
- Graves HC, Wehner N, Stamey TA. Comparison of a polyclonal and monoclonal immunoassay for PSA: need for an international antigen standard. *J Urol*. 1990;144 (6) : 1516-1522.
- Lilja H, Christensson A, Dahlen U, et al. Prostate-specific antigen in human serum occurs predominantly in complex with alpha 1-antichymotrypsin. *Clin Chem*. 1991;37 (9) :1618-1625.
- Gormley GJ, Stoner E, Bruskewitz RC, et al. The effect of finasteride in men with benign prostatic hyperplasia. The Finasteride Study Group. *N Engl J Med*. 1992;327 (17) :1185-1191.
- Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45 (7) :942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F (ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38 (9) :1737-1742.
- Wang MC, Valenzuela LA, Murphy GP, Chu TM. Purification of a human prostate specific antigen. *Invest Urol*. 1979;17 (2) :159-163.
- Watt KW, Lee PJ, M'Timkulu T, et al. Human prostate-specific antigen: structural and functional similarity with serine proteases. *Proc Natl Acad Sci USA*. 1986;83 (10) :3166-3170.
- Lilja H. A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. *J Clin Invest*. 1985;76 (5) :1899-1903.
- Stamey TA, Yang N, Hay AR, et al. Prostate-specific antigen as a serum marker for adenocarcinoma of the prostate. *N Engl J Med*. 1987;317 (15) :908-916.
- Schellhammer PF, Schlossberg SM, el-Mahdi AM, et al. Prostate specific antigen levels after definitive irradiation for carcinoma of the prostate. *J Urol*. 1991;145 (5) :1008-1010.

- Cooner WH, Mosley BR, Rutherford CL Jr, et al. Prostate cancer detection in a clinical urological practice by ultrasonography, digital rectal examination and prostate specific antigen. 1990. *J Urol*. 2002;167 (2, pt 2) :966-973.
- Oesterling JE. Prostate specific antigen: a critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. *J Urol*. 1991;145 (5) :907-923.
- Stamey TA, Kabalin JN, McNeal JE, et al. Prostate specific antigen in the diagnosis and treatment of adenocarcinoma of the prostate. II. Radical prostatectomy treated patients. *J Urol*. 1989;141 (5) :1076-1083.
- Lange PH, Ercole CJ, Lightner DJ, et al. The value of serum prostate specific antigen determinations before and after radical prostatectomy. *J Urol*. 1989;141 (4) :873-879.
- Killian CS, Yang N, Emrich LJ, et al. Prognostic importance of prostate-specific antigen for monitoring patients with stages B2 to D1 prostate cancer. *Cancer Res*. 1985;45 (2) :886-891.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- **26. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.
- Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
- Clinical and Laboratory Standards Institute (formerly NCCLS). *Primary reference preparations used to standardize calibration of immunochemical assays for serum prostate specific antigen (PSA) ; Approved Guideline*. Wayne, PA: Clinical and Laboratory Standards Institute; 1997. NCCLS Document I/LA19-A.

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