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- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition*. NCCLS document H3-A5 [ISBN 1-56238-515-1]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Kaplan, L.A, Pesce, A.J.: *Clinical Chemistry Theory Analysis and Correlation*, third edition, Mosby Co, St. Louis, MO, 1996: p 552.
- Burtis, C.A, Ashwood, E.R: *Tietz Textbook of Clinical Chemistry*, WB Saunders Co, Philadelphia, PA, 1999, p 42-56, 1407 (Specimen Collection and Handling of Specimens for Testing), p.1829 (Pathophysiology).
- Clinical and Laboratory Standards Institute /NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition*. NCCLS document H1-A5 [ISBN 1-56238-519-4]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Fiske, CH, Subbarow, U, The colorimetric determination of phosphorus, *J Biol Chem* 1925; 66:375.
- Gomori, G, A modification of the colorimetric phosphorus determination for use with the photoelectric colorimeter, *J Lab Clin Med* 1941; 27:955.
- Drewes, P, Direct colorimetric determination of phosphorus in serum and urine, *Clin Chem Acta* 1972; 39:81.
- Goodwin, JF, Quantification of serum inorganic phosphorus, phosphatase, and urinary phosphate without preliminary treatment, *Clin Chem* 1970; 16:776.
- Clinical and Laboratory Standards Institute /NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests*, 3rd ed. Washington, DC: AACC Press, 2007, p 720.

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