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- Ullman EF, Kirakossian H, Switchenko AC, Ishkanian J, et. al., Luminescent oxygen channeling assay (LOCI™): sensitive, broadly applicable homogenous immunoassay method. Clin Chem 42:9 1996, 1518-1526.
- Ullman EF, Kirakossian H, Sharat S, Ping Wu Z, Irvin BR, et. al Luminescent oxygen channeling immunoassay: Measurement of particle binding kinetics by chemiluminescence. Proc.Natl.Acad.Sci. USA, Vol 91, pp.5426-5430, June 1994 Biochemistry.
- Tietz NW, Textbook of Clinical Chemistry, W.B. Saunders Co. Philadelphia, PA 1999; pp. 42-72.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Jung K, et al. Preanalytical Determinants of Total and Free Prostate-Specific Antigen and Their Ratio: Blood Collection and Storage Conditions. Clin Chem 1998; 44: 685-688.
- Crawford ED, Schutz MJ, Clejan S, et al. The effect of digital rectal examination of prostate specific antigen levels. JAMA 1992; 267: 2227.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Fraser, Callum G. PhD, Biological Variation: From Principles to Practice, Washington, DC: AACCPress, 2001 (pp. 71-78).
- Ricos C, Alvarez V, Cava F, Garcia-Lario JV, Hernandez A, Jimenez CV, Minchinela J, Perich C, Simon M. Current databases on biologic variation: pros, cons, and progress. Scand J Clin Lab Invest, 1999; 59: 491-500.
- Ryall RG, Story CJ, and Turner DR. Reappraisal of the causes of the “hook”effect in two-site immunometric assays. Anal Biochem 1982; 127: 308-15.
- Morgan WR, Zincke H, et al. Prostate specific antigen values after radical retropubic prostatectomy for adenocarcinoma of the prostate: Impact of adjuvant treatment (hormonal and radiation).J Urol 1991; 145: 319.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F(ab')2 conjugate and polyclonal mouse IgG. Clin Chem 1992; 38: 1737-42.
- Kricka LJ. Human Anti-Animal Antibody Interferences in Immunological Assays. Clin Chem 1999; 45:7; 942-956.
- Lilja H. A kallikrein-like serine protease in prostatic fluid cleaves the predominant seminal vesicle protein. J Clin Invest 1985;76: 1899-903.
- Watt KWK, Lee P-J, Timkulu T, Chan W-P, Loor R. Human prostate-specific antigen: structural and functional similarities with serine protease. Proc Natl Acad Sci. 1986; 83: 3166-3170.
- Coley CM, Barry MJ, Fleming C, et al: Early detection of prostate cancer: I. Prior probability and effectiveness of tests. ANN Intern Med 1997; 126: 394-406.
- Stenman V-H, Leinonen J, Althtan H, et al. A complex between prostate-specific antigen and α_1 -antichymotrypsin is the major form of prostate-specific antigen in serum of patients with prostate cancer: Assay of the complex improves clinical sensitivity for cancer. Cancer Res 1991;51: 222-6.
- Lilja H, Christensson A, Dahlen V, et al. Prostate-specific antigen in human serum occurs predominantly in complex with α_1 -antichymotrypsin. Clin Chem 1991; 37: 1618-25.
- www.cancer.org
- Rogers E, Ohori M, Kassabian VS, et al: Salvage radical prostatectomy: Outcome measured by serum prostate specific antigen levels. J Urol 1995; 153: 104-110.
- Pound CR, Partin AW, Eisenberger MA, et al: Natural history of progression after PSA elevation following radical prostatectomy. JAMA 1999; 281: 1591-1597.
- Morgan WR, Zincke H, et al. Prostate specific antigen values after radical retropubic prostatectomy for adenocarcinoma of the prostate: Impact of adjuvant treatment (hormonal and radiation). J Urol 1991; 145: 319.
- Lee W, Hanlon AL, Hanks GE: Prostate specific antigen nadir level and disease-free survival. J.Urol 1996; 156: 450-453.
- Clinical and Laboratory Standards Institute/NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Roddam AW, Rimmer J, Nickerson C, Ward AM. Prostate-specific antigen: bias and molarity of commercial assays for PSA in use in England. Ann Clin Biochem 2006;43:35-48.

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