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1. Shaw LM, Nicholls A, Hale M, et al. Therapeutic monitoring of mycophenolic acid: A consensus panel report. *Clin Biochem*. 1998;31/5:317-322.
2. Nicholls AJ. Opportunities for therapeutic monitoring of mycophenolate mofetil dose in renal transplantation suggested by the pharmacokinetic/pharmacodynamic relationship for mycophenolic acid and suppression of rejection. *Clin Biochem*. 1998;31/5:329-333.
3. Tsina I, Chu F, Hama K, et al. Manual and automated (robotic) high-performance liquid chromatography methods for the determination of mycophenolic acid and its glucuronide conjugate in human plasma. *J of Chromatogr B*. 1996;675:119-129.
4. JY Richmond, RW McKinney, eds. *Biosafety in microbiological and biomedical laboratories*. 3rd ed. Washington, DC: Centers for Disease Control and Prevention; 1993. US Dept of Health and Human Services publication 93-8395.
5. Bloodborne Pathogens (29 CFR; Ch XVII:Part 1910-1030) *Federal Register*. July 1, 1991.
6. Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
7. European Mycophenolate Mofetil Cooperative Study Group. Placebo-controlled study of mycophenolate mofetil combined with cyclosporin and corticosteroids for prevention of acute rejection. *Lancet*. 1995;345/8961:1321-1325.
8. Sollinger HW for the US Renal Transplant Mycophenolate Mofetil Study Group. Mycophenolate mofetil for the prevention of acute rejection in primary cadaveric renal allograft recipients. *Transplantation*. 1995;60/3:225-232.
9. Kobashigawa J, Miller L, Renlund D, et al. A randomized active-controlled trial of mycophenolate mofetil in heart transplant recipients. *Transplantation*. 1998;66/4:507-515.
10. Allison AC, Kowalski WJ, Muller CD, et al. Mechanisms of action of mycophenolic acid. In: Allison AC, Lafferty K, Fliri H, eds. *Immunosuppressive and Antiinflammatory Drugs*. New York, NY: The New York Academy of Sciences; 1993;696:63-87.
11. Wu JC. Mycophenolate mofetil: Molecular mechanisms of action. *Perspect Drug Discovery Design*. 1994;2:185-204.
12. Bullingham RES, Nicholls AJ, Kamm BR. Clinical pharmacokinetics of mycophenolic mofetil. *Clin Pharmacokinet*. 1998;34/6:429-455.
13. Schütz E, Shipkova M, Armstrong VW, et al. Therapeutic drug monitoring of mycophenolic acid: comparison of HPLC and immunoassay reveals new MPA metabolites. *Transplant Proc*. 1998;30:1185-1187.
14. Shipkova M, Armstrong VW, Wieland E, Niedmann PD, Schütz E, Brenner-Weiss G, Voihsel M, Braun F, Oellerich M. Identification of glucoside and carboxyl-linked glucuronide conjugates of mycophenolic acid in plasma of transplant recipients treated with mycophenolate mofetil. *Br J Pharmacol*. 1999;126/5:1075-1082.
15. Schütz E, Shipkova M, Armstrong VW, Wieland E, Oellerich M. Identification of a pharmacologically active metabolite of mycophenolic acid in plasma of transplant recipients treated with mycophenolate mofetil. *Clin Chem*. 1999;45/3:419-422.
16. Nowak I, Shaw LM. Mycophenolic acid binding to human serum albumin: Characterization and relation to pharmacodynamics. *Clin Chem*. 1995;41/7:1011-1017.
17. Shaw LM, Nowak I. Mycophenolic acid: measurement and relationship to pharmacologic effects. *Ther Drug Monit*. 1995;17/6:685-689.
18. Meier-Kriesche H-U, Shaw LM, Korecka M, Kaplan B. Pharmacokinetics of mycophenolic acid in renal insufficiency. *Ther Drug Monit*. 2000;22/1:27-30.
19. Weber LT, Shipkova M, Lamersdorf T, et al. Pharmacokinetics of mycophenolic acid (MPA) and determinants of MPA free fraction in pediatric and adult renal transplant recipients. *J of Am Soc Nephrol*. 1998;9/8:1511-1520.
20. Weber LT, Lamersdorf T, Shipkova M, Niedmann PD, Wiesel M, Zimmerhackl LB, Staskewitz A, Schütz E, Mehls O, Oellerich M, Armstrong VW, Tönshoff B, members of the German study group on mycophenolate mofetil therapy in pediatric renal transplant recipients. Area under the plasma concentration-time curve for total, but not free, mycophenolic acid increases in the stable phase after renal transplantation: a longitudinal study in pediatric patients. *Ther Drug Monit*. 1999;21:498-506.
21. Oellerich M, Shipkova M, Schütz E, Wieland E, Weber L, Tönshoff B, Armstrong VW, members of the German study group on mycophenolate mofetil therapy in pediatric renal transplant recipients. Pharmacokinetic and metabolic investigations of mycophenolic acid in pediatric patients after renal transplantation: implications for therapeutic drug monitoring. *Ther Drug Monit*. 2000;22/1:20-26.
22. Weber LT, Schütz E, Lamersdorf T, Shipkova M, Niedmann PD, Oellerich M, Zimmerhackl LB, Staskewitz A, Mehls O, Armstrong VW, Tönshoff B, German study group on mycophenolate mofetil (MMF) therapy. Therapeutic drug monitoring of total and free mycophenolic acid (MPA) and limited sampling strategy for determination of MPA-AUC in paediatric renal transplant recipients. *Nephrol Dial Transplant*. 1999;14/Suppl4:34-35.
23. Shaw LM, Kaplan B, DeNofrio D, Korecka M, Brayman KL. Pharmacokinetics and concentration-control investigations of mycophenolic acid in adults after transplantation. *Ther Drug Monit*. 2000;22/1:14-19.
24. RS-61443 Investigation Committee. Pilot study of mycophenolate mofetil (RS-61443) in the prevention of acute rejection following renal transplantation in Japanese patients. *Transplant Proc*. 1995;27:1421-1424.
25. van Gelder T, Hilbrands LB, Vanrenterghem Y, et al. A randomized double-blind, multicenter plasma concentration controlled study of the safety and efficacy of oral mycophenolate mofetil for the prevention of acute rejection after kidney transplantation. *Transplantation*. 1999;68/2:261-266.
26. Hale MD, Nicholls AJ, Bullingham RES, Hené R, Hoitsma A, Squifflet J-P, Weimar W, Vanrenterghem Y, Van de Woude FJ, Verpooten GA. The pharmacokinetic-pharmacodynamic relationship for mycophenolate mofetil in renal transplantation. *Clin Pharmacol Ther*. 1998;64:672-683.
27. Shaw LM, Korecka M, Aradhye S, et al. Analysis of mycophenolic acid pharmacokinetics in African American and Caucasian renal transplant recipients. *Transplantation*. 1998;66/8:S5.
28. Zanker B, Sohr B, Eder M, Frohmann E, Land W. Comparison of MPA trough levels in patients with severe diabetes mellitus and from non-diabetics after transplantation. *Transplant Proc*. 1999;31:1167.
29. Morgera S, Budde K, Lampe D, Ahnert V, Fritsche L, Kuchinke S, Neumayer H-H. Mycophenolate mofetil pharmacokinetics in renal transplant recipients on peritoneal dialysis. *Transpl Int*. 1998;11:53-57.
30. Wollenberg K, Krumme B, Pisarski P, Schollmeyer P, Kirste G. Pharmacokinetics of mycophenolic acid in the early period after kidney transplantation. *Transplant Proc*. 1998;30:4090-4091.
31. Jacqz-Aigrain E, Shaghghi K, Baudouin V, Popon M, Zhang D, Maisin A, Loirat C. Pharmacokinetics and tolerance of mycophenolate mofetil in renal transplant children. *Pediatr Nephrol*. 2000;14/2:95-99.
32. Schütz E, Andreeva M, Niedmann PD, et al. Mycophenolic acid determination after solid organ transplantation: RP-HPLC and Emit® compared. *Ther Drug Monit*. 1997;19/5:566.
33. Dias VC, Martens C, Zimmermann L, et al. Evaluation of the Emit® Mycophenolic Acid Assay on the COBAS MIRA®.

- Ther Drug Monit.* 1997;19/5:563.
34. Vogl M, Weigel G, Seebacher G, et al. Evaluation of the Emit® mycophenolic acid assay from Siemens. *Ther Drug Monit.* 1999;21/6:638-643.
35. Beal JL, Jones CE, Taylor PJ, Tett SE. Evaluation of an immunoassay (Emit®) for mycophenolic acid in plasma from renal transplant recipients compared with a high-performance liquid chromatography assay. *Ther Drug Monit.* 1999;20:685-690.
36. Yeung JS, Wang W, Chan L. Determination of mycophenolic acid level: comparison of high-performance liquid chromatography with homogenous enzyme-immunoassay. *Transplant Proc.* 1999;31:1214-1215.
37. Brunet M, Oppenheimer F, Martorell J, Vilardell J, Carreño MC, Carrillo M, Corbella J. Mycophenolic acid monitoring: evaluation of the Emit® MPA immunoassay in kidney and lung transplantation. *Transplant Proc.* 1999;31:2275-2276.
38. Shipkova M, Schütz E, Armstrong VW, Niedmann PD, Wieland E, Oellerich M. Overestimation of mycophenolic acid by Emit® correlates with MPA metabolite. *Transplant Proc.* 1999;31:1135-1137.
39. Johnson AG, Rigby RJ, Taylor PJ, Jones CE, Allen J, Franzen K, Falk MC, Nicol D. The kinetics of mycophenolic acid and its glucuronide metabolite in adult kidney transplant recipients. *Clin Pharmacol Ther.* 1999;66:492-500.
40. Krumme B, Wollenberg K, Kirste G, Schollmeyer P. Drug monitoring of mycophenolic acid in the early period after renal transplantation. *Transplant Proc.* 1998;30:1773-1774.
41. Meiser BM, Pfeiffer M, Schmidt D, et al. The efficacy of the combination of tacrolimus and mycophenolate mofetil for prevention of acute myocardial rejection is dependent on routine monitoring of mycophenolic acid trough acid levels. *Transplant Proc.* 1999;31:84-87.
42. Schütz E, Armstrong VW, Shipkova M, et al. Limited sampling strategy for the determination of mycophenolic acid area under the curve in pediatric kidney recipients. *Transplant Proc.* 1998;30:1182-1184.
43. Sanquer S, Breil M, Baron C, et al. Induction of inosine monophosphate dehydrogenase activity after long-term treatment with mycophenolate mofetil. *Clin Pharmacol Ther.* 1999;65:640-648.
44. Zanker B, Schneeberger H, Rothenpieler U, et al. Mycophenolate mofetil-based, cyclosporine-free induction and maintenance immunosuppression. *Transplantation.* 1998;66/1:44-49.
45. Smak Gregoor PJH, de Sévaux RGL, Hené RJ, Hesse CJ, Hilbrands LB, Vos P, van Gelder T, Hoitsma AJ, Weimar W. Effect of cyclosporine on mycophenolic acid trough levels in kidney transplant recipients. *Transplantation.* 1999;68/10:1603-1606.
46. Zucker K, Tsaroucha A, Olson L, Esquenazi V, Tzakis A, Miller J. Evidence that tacrolimus augments the bioavailability of mycophenolate mofetil through the inhibition of mycophenolic acid glucuronidation. *Ther Drug Monit.* 1999;21:35-43.
47. Hübner G, Eismann R, Sziegoleit W. Drug interaction between mycophenolate mofetil and tacrolimus detectable within therapeutic mycophenolic acid monitoring in renal transplant patients. *Ther Drug Monit.* 1999;21:536-539.
48. Filler G, Zimmering M, Mai I. Pharmacokinetics of mycophenolate mofetil are influenced by concomitant immunosuppression. *Pediatr Nephrol.* 2000;14/2:100-104.
49. Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition.* CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
50. *Clinical and Laboratory Standards Institute (formerly*

NCCLS). Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline—Second Edition. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

51. Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline.* CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.

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