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- Tokuda K, Tanimoto K. New method of measuring serum bilirubin using vanadic acid. *Jpn J Clin Chem.* 1993;22 (2) :116-122.
- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Burtis CA, Ashwood ER, Bruns DE, Tietz NW, *Textbook of Clinical Chemistry and Molecular Diagnostics.* 4th ed., Elsevier Saunders, 2006:1165-1208.
- Young DS, *Effects of Preanalytical Variables on Clinical Laboratory Tests.* 2nd ed. Washington, DC: AACC Press, 1997:3-88,4-492.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Data on file at Siemens Healthcare Diagnostics.
- Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Tietz NW. *Clinical Guide to Laboratory Tests.* 3rd ed. Philadelphia, PA: WB Saunders Company; 1995: 88-91.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests.* 3rd ed. Washington, DC: AACC Press; 1990.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c-ed3.
- Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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