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- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
- Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45 (7) :942-956.
- Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F (ab')₂ conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38 (9) :1737-1742.
- Austin L, Heath H. Calcitonin: physiology and pathophysiology. *N Engl J Med*. 1981; 304 (5) :269-278.
- Zaidi M, Moonga BS, Bevis PJ, et al. Expression and function of the calcitonin gene products. *Vitam Horm*. 1991;46:87-164.
- Faour O, Gilloteaux J. Calcitonin: Survey of new anatomy data to pathology and therapeutic aspects. *Translational Research in Anatomy*. 2017;6:4-15.
- Inzerillo AM, Ziadi M, Huang CL. Calcitonin: physiological actions and clinical applications. *J Pediatr Endocrinol Metab*. 2004; 17 (7) :931-940.
- Boron WF, Boulpaep EL. *Medical Physiology: A Cellular And Molecular Approach*. Philadelphia, PA: Saunders; 2004.
- Guesgen C, Willms A, Zwad A, et al. Investigation of factors potentially influencing calcitonin levels in the screening and follow-up for medullary thyroid carcinoma: a cautionary note. *BMC Clinical Pathology*. 2013; Nov 4; 13 (1) :27:1-8.
- Society for Endocrinology. You and your hormones; calcitonin. 2011. <http://www.yourhormones.info/hormones/calcitonin.aspx>. Reviewed Jan 2015. Accessed June 20, 2017.
- Hahn JR, Lee MS, Min YK, et al. Routine measurement of serum calcitonin is useful for early detection of medullary thyroid carcinoma in patients with nodular thyroid diseases. *Thyroid*. 2001; 11 (1) :73-80.

- Pacini F, Fontanelli M, Fugazzola L, et al. Routine measurement of serum calcitonin in nodular thyroid diseases allows the preoperative diagnosis of unsuspected sporadic medullary thyroid carcinoma. *J Clin Endocrinol Metab*. 1994 Apr;78 (4) :826-829.
- Hillyard C, Oscier DG, Foa R, et al. Immunoreactive calcitonin in leukaemia. *Br Med J*. 1979 Dec 1;2 (6202) :1392-1393.
- Hillyard C, Stevenson J. Calcitonin as a tumour marker. *Ligand Rev*. 1980 (Spring) ;2 (1) :31-33.
- Milhaud G. Ectopic secretion of calcitonin. In: *Calcitonin 1980; Chemistry, Physiology, Pharmacology, and Clinical Aspects*. Pecile A, ed. Amsterdam: *Excerpta Medica*. 1981: 154-169.
- Mulder H, Hackeng WHL. Ectopic secretion of calcitonin. *Acta Med Scand*. 1978; 204: 253-256.
- Schwartz KE, Wolfsen AR, Forster B, et al. Calcitonin in nonthyroidal cancer. *J Clin Endocrinol Metab*. 1979 Sep;49 (3) :438-444.
- Wallach S, Royston I, Taetle R, et al. Plasma calcitonin as a marker of disease activity in patients with small cell carcinoma of the lung. *J Clin Endocrinol Metab*. 1981; Sep;53 (3) :602-606.
- Mulder H, Silberbusch J, Hackeng WH, et al. Hypercalcitoninaemia in patients with chronic inflammatory disease. *Neth J Med*. 1980;23 (3) :129-131.
- Toffaletti J. Calcitonin. *Clin Chem News*. 1984 (July) ;10 (7) :20-21.
- Williams ED, Brown CL, Doniach I. Pathological and clinical findings in a series of 67 cases of medullary carcinoma of the thyroid. *J Clin Pathol*. 1966; Mar;19 (2) :103-113.
- Hadoux J, Schlumberger M. Chemotherapy and tyrosine-kinase inhibitors for medullary thyroid cancer. *Best Pract Res Clin Endocrinol Metab*. 2017 Jun;31 (3) :335-347.
- Leboulleux S, Baudin E, Travagli JP, et al. Medullary thyroid carcinoma. *Clin Endocrinol (Oxf)*. 2004 Sep;61 (3) :229-310.
- Elisei R, Pinchera A. Advances in the follow-up of differentiated or medullary thyroid cancer. *Nat Rev Endocrinol*. 2012;8:466-475.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
- Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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