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2×10.0 mL

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1. Cristancho E, Riveros A, Sanchez A, et al. Diurnal changes of arterial oxygen saturation and erythropoietin concentration in male and female highlanders. *Physiol Rep*. Sep;4(17).
2. Wide L, Bengtsson C, Birgegård G. Circadian rhythm of erythropoietin in human serum. *Br J Haematol*. 1989 May;72(1):85-90.
3. Cahan C, Decker MJ, Arnold JL, et al. Diurnal variations in serum erythropoietin levels in healthy subjects and sleep apnea patients. *J Appl Physiol* (1985). 1992 Jun;72(6):2112-2117.
4. Miller M, et al. Diurnal levels of immunoreactive erythropoietin in normal subjects and subjects with lung disease. *British J of Haematology*. 1981, 49, 189-200.
5. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard - Sixth Edition*. CLSI Document GP41-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
6. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard- Sixth Edition*. CLSI Document GP39-A6. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
7. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline - Fourth Edition*. CLSI Document GP44-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
8. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guidelines - Second Edition*. CLSI Document EP07-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
9. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline - Third Edition*. CLSI Document EP28-A3c. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
10. Boscato LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. *Clin Chem*. 1988;34:27-33.
11. Naderi N, Yang DT. Lymphoplasmacytic lymphoma and Waldenström macroglobulinemia. *Arch Pathol Lab Med*. 2013 Apr;137(4):580-585.
12. Eckardt KU, Bauer C. Erythropoietin in health and disease. *Europ J Clin Invest*. 1989 Apr;19(2):117-127.
13. Jelkmann W. Erythropoietin: structure, control of production, and function. *Physiol Rev*. 1992 Apr;72(2):449-489.
14. Jelkmann W. Regulation of erythropoietin production. *J Physiol*. 2011 Mar 15;589(Pt 6):1251-8.
15. Tsiolakidou G, Koutroubakis IE. Stimulating erythropoiesis in inflammatory bowel disease associated anemia. *World J Gastroenterol*. 2007 Sep;13(36):4798-4806.
16. Pulte ED, McKenzie SE, Caro J, Ballas SK. Erythropoietin levels in patients with sickle cell disease do not correlate with known inducers of erythropoietin. *Hemoglobin*. 2014;38(6):385-389.
17. Gaya A, Urbano-Ispizua A, Fernández-Avilés F, Salamero O, et al. Anemia associated with impaired erythropoietin secretion after allogeneic stem cell transplantation: incidence, risk factors, and response to treatment. *Biol Blood Marrow Transplant*. 2008 Aug;14(8):880-887.
18. Recny MA, Scoble HA, Kim Y. Structural characterization of natural human urinary and recombinant DNA-derived erythropoietin. Identification of des-arginine 166 erythropoietin. *J Biol Chem*. 1987 Dec 15;262(35):17156-17163.
19. Dordal MS, Wang FF, Goldwasser E. The role of carbohydrate in erythropoietin action. *Endocrinology*. 1985 Jun;116(6):2293-2299.
20. Miyake T, Kung CK-H, Goldwasser. Purification of human erythropoietin. *J Bio Chem*. 252,(15), 5558-5564. 1977.
21. Bunn HF. Erythropoietin. *Cold Spring Harb Perspect Med*. 2013 Mar 1;3(3):a011619.
22. Eschbach JW, Adamson JW. Recombinant human erythropoietin: implications for nephrology. *Am J Kidney Dis*. 1988; 11:203-209.
23. Koch KM, Kuhn K, Nonnast-Daniel B, Scigalla P, eds. *Treatment of renal anemia with recombinant human erythropoietin*. In: Berlyne GM, Giovannetti S, series editors. *Contributions to Nephrology*. New York: Karger, 1988; 66:1-15 and 54-61.

24. Gross AW, Lodish HF. Cellular trafficking and degradation of erythropoietin and novel erythropoiesis stimulating protein (NESP). *J Biol Chem.* 2006 Jan 27;281(4):2024-2032.
25. Eschbach JW, Kelly MR, Haley NR, et al. Treatment of the anemia of progressive renal failure with recombinant human erythropoietin. *N Engl J Med.* 1989 Jul 20;321(3):158-163.
26. Haase VH. Hypoxic regulation of erythropoiesis and iron metabolism. *Am J Physiol Renal Physiol.* 2010 Jul;299(1):F1-13.
27. Egrie JC, Strickland TW, Lane J, et al. Characterization and biological effects of recombinant human erythropoietin. *Immunobiology.* 1986 Sep;172(3-5):213-224.
28. Lin FK, et al. Cloning and expression of the human erythropoietin gene. *Proc Natl Acad Sci U S A.* 1985 Nov;82(22):7580-7584.
29. Egrie JC, Browne JK. Development and characterization of novel erythropoiesis stimulating protein (NESP). *Br J Cancer.* 2001 Apr;84 Suppl 1:3-10.
30. Owen WE, Roberts W. Performance characteristics of a new Immulite 2000 system erythropoietin assay. *Clin Chim Acta.* 2011 Feb 20;412 (5-6):480-482.
31. Messinezy M, Westwood NB, El-Hemaidi I, et al. Serum erythropoietin values in erythrocytoses and in primary thrombocythaemia. *Br J Haematol.* 2002 Apr;117(1):47-53.
32. Mossuz P, Girodon F, Hermouet S, et al. Serum erythropoietin measured by chemiluminescent immunometric assay: an accurate diagnostic test for absolute erythrocytosis. *Clin Chem.* 2005 Jun;51 (6):1018-1021.
33. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline - Third Edition.* CLSI Document EP09-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2013.
34. Clinical and Laboratory Standards Institute. *Evaluation of the Linearity of Quantitative Measurement Procedures. A Statistical Approach; Approved Guideline.* CLSI Document EP06-A. Wayne, PA: Clinical and Laboratory Standards Institute; 2003.
35. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guidelines - Third Edition.* CLSI Document EP05-A3. Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
36. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline - Second Edition.* CLSI Document EP17-A2. Wayne, PA: Clinical and Laboratory Standards Institute; 2012.

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