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1. Cristancho E, Riveros A, Sanchez A, et al. Diurnal changes of arterial oxygen saturation and erythropoietin concentration in male and female highlanders. *Physiol Rep*. Sep;4 (17).
2. Wide L, Bengtsson C, Birgegard G. Circadian rhythm of erythropoietin in human serum. *Br J Haematol*. 1989 May;72 (1) :85-90.
3. Cahan C, Decker MJ, Arnold JL, et al. Diurnal variations in serum erythropoietin levels in healthy subjects and sleep apnea patients. *J Appl Physiol* (1985). 1992 Jun;72 (6) :2112-2117.
4. Miller M, et al. Diurnal levels of immunoreactive erythropoietin in normal subjects and subjects with lung disease. *British J of Haematology*. 1981, 49, 189-200.
5. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers from Occupationally Acquired Infections; Approved Guideline-Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
6. Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
7. Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard-Sixth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
8. Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline-Fourth Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
9. Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
10. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
11. Naderi N, Yang DT. Lymphoplasmacytic lymphoma and Waldenström macroglobulinemia. *Arch Pathol Lab Med*. 2013 Apr;137 (4) :580-585.
12. Kricka LJ. Human anti-animal antibody interferences in immunological assays. *Clin Chem*. 1999;45 (7) :942-956.

13. Vaidya HC, Beatty BG. Eliminating interference from heterophilic antibodies in a two-site immunoassay for creatine kinase MB by using F (ab')<sub>2</sub> conjugate and polyclonal mouse IgG. *Clin Chem*. 1992;38 (9) : 1737-1742.
14. Eckardt KU, Bauer C. Erythropoietin in health and disease. *Europ J Clin Invest*. 1989 Apr;19 (2) :117-127.
15. Jelkmann W. Erythropoietin: structure, control of production, and function. *Physiol Rev*. 1992 Apr;72 (2) :449-489.
16. Jelkmann W. Regulation of erythropoietin production. *J Physiol*. 2011 Mar 15;589 (Pt 6) :1251-8.
17. siolakidou G, Koutroubakis IE. Stimulating erythropoiesis in inflammatory bowel disease associated anemia. *World J Gastroenterol*. 2007 Sep;13 (36) :4798-4806.
18. Pulte ED, McKenzie SE, Caro J, Ballas SK. Erythropoietin levels in patients with sickle cell disease do not correlate with known inducers of erythropoietin. *Hemoglobin*. 2014;38 (6) :385-389.
19. Gaya A, Urbano-Ispizua A, Fernández-Avilés F, Salamero O, et al. Anemia associated with impaired erythropoietin secretion after allogeneic stem cell transplantation: incidence, risk factors, and response to treatment. *Biol Blood Marrow Transplant*. 2008 Aug;14 (8) :880-887.
20. Recny MA, Scoble HA, Kim Y. Structural characterization of natural human urinary and recombinant DNA-derived erythropoietin. Identification of des-arginine 166 erythropoietin. *J Biol Chem*. 1987 Dec 15;262 (35) :17156-17163.
21. Dordal MS, Wang FF, Goldwasser E. The role of carbohydrate in erythropoietin action. *Endocrinology*. 1985 Jun;116 (6) :2293-2299.
22. Miyake T, Kung CK-H, Goldwasser. Purification of human erythropoietin. *J Bio Chem*. 252, (15), 5558-5564. 1977.
23. Bunn HF. Erythropoietin. *Cold Spring Harb Perspect Med*. 2013 Mar 1;3 (3) :a011619.
24. Eschbach JW, Adamson JW. Recombinant human erythropoietin: implications for nephrology. *Am J Kidney Dis*. 1988; 11:203-209.
25. Koch KM, Kuhn K, Nonnast-Daniel B, Scigalla P, eds. *Treatment of renal anemia with recombinant human erythropoietin*. In: Berlyne GM, Giovannetti S, series editors. Contributions to Nephrology. New York: Karger, 1988; 66:1-15 and 54-61.
26. Gross AW, Lodish HF. Cellular trafficking and degradation of erythropoietin and novel erythropoiesis stimulating protein (NESP). *J Biol Chem*. 2006 Jan 27;281 (4) :2024-2032.
27. Eschbach JW, Kelly MR, Haley NR, et al. Treatment of the anemia of progressive renal failure with recombinant human erythropoietin. *N Engl J Med*. 1989 Jul 20;321 (3) :158-163.
28. Haase VH. Hypoxic regulation of erythropoiesis and iron metabolism. *Am J Physiol Renal Physiol*. 2010 Jul;299 (1) :F1-13.
29. Egrie JC, Strickland TW, Lane J, et al. Characterization and biological effects of recombinant human erythropoietin. *Immunobiology*. 1986 Sep;172 (3-5) :213-224.
30. Lin FK, et al. Cloning and expression of the human erythropoietin gene. *Proc Natl Acad Sci U S A*. 1985 Nov;82 (22) :7580-7584.
31. Egrie JC, Browne JK. Development and characterization of novel erythropoiesis stimulating protein (NESP). *Br J Cancer*. 2001 Apr;84 Suppl 1:3-10.
32. Owen WE, Roberts W. Performance characteristics of a new Immulite 2000 system erythropoietin assay. *Clin Chim Acta*. 2011 Feb 20;412 (5-6) :480-482.
33. Messinezy M, Westwood NB, El-Hemaidi I, et al. Serum erythropoietin values in erythrocytosis and in primary thrombocythaemia. *Br J Haematol*. 2002 Apr;117 (1) :47-53.
34. Mossuz P, Girodon F, Hermouet S, et. al. Serum erythropoietin measured by chemiluminescent immunometric assay: an accurate diagnostic test for absolute erythrocytosis. *Clin Chem*. 2005 Jun;51 (6) :1018-1021.
35. Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
36. Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
37. Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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