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1. Kamal AH, Tefferi A, Pruthi RK. How to interpret and pursue an abnormal prothrombin time, activated partial thromboplastin time, and bleeding time in adults. Mayo Clin Proc 2007; 82: 864-73.
2. Leung LLK. Perioperative evaluation of bleeding diathesis. Hematology Am Soc Hematol Educ Program 2006; 457-61.

3. Eikelboom JW, Hirsh J. Monitoring unfractionated heparin with the aPTT: Time for a fresh look. *Thromb Haemost* 2006; 96: 547-52

4. Bates SM, Johnston M, Hirsh J, Ginsberg JS. Use of a fixed activated partial thromboplastin time ratio to establish a therapeutic range for unfractionated heparin. *Arch Intern Med* 2001; 161: 385-91.

5. Hirsh J, Raschke R. Heparin and low-molecular-weight heparin. The seventh ACCP Conference on Antithrombotic and Thrombolytic Therapy. *Chest* 2004; 126: 188S-203S.

6. Olson JD, Arkin CF, Brandt JT, et al. College of American Pathologists Conference XXXI on laboratory monitoring of anticoagulant therapy: laboratory monitoring of unfractionated heparin therapy. *Arch Pathol Lab Med.* 1998; 122: 782-98

7. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline—Fifth Edition. CLSI document H21-A5 (ISBN 1-56238-657-3). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2008.

8. Hattersley PG, Hayse D. The effect of increased contact activation time on the activated partial thromboplastin time. *Am J Clin Pathol* 1976; 66: 479-82.

9. Ambrus JL, Schimert G, Lajos TZ, et al. Effect of antifibrinolytic agents and estrogens on blood loss and blood coagulation factors during open heart surgery. *J Med* 1971; 2: 65-81.

10. Prasad RN, Koh SC, Viegas OA, Ratnam SS. Effects on hemostasis after two-year use of low dose combined oral contraceptives with gestodene or levonorgestrel. *Clin Appl Thromb Hemost* 1999; 51: 60-70.

11. Solomon GE, Hilgartner MW, Kutt H. Coagulation defects caused by diphenylhydantoin. *Neurology* 1972; 22: 1165-71.

12. Young DS, ed. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington: AACC Press, 2000.

13. Greinacher A, Warkentin TE. The direct thrombin inhibitor hirudin. *Thromb Haemost* 2008; 99: 819-29.

14. Soloway HB, Cox SP, Donahoo JV. Sensitivity of the activated partial thromboplastin time to heparin: effect of anticoagulant used for sample collection. *Am J Clin Pathol* 1973; 59: 760-2.

15. Lippi G, Franchini M, Montagnana M et al. Quality and reliability of routine coagulation testing: can we trust that sample? *Blood Coag Fibrinol* 2006; 17: 513-9.

16. CLSI. One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline. CLSI document H47-A [ISBN 1-56238-301-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 1996.

17. CLSI. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline - Second Edition. CLSI document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2000.

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