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1. Kamal AH, Tefferi A, Pruthi RK. How to interpret and pursue an abnormal prothrombin time, activated partial thromboplastin time, and bleeding time in adults. Mayo Clin Proc 2007; 82: 864-73.
2. Leung LLK. Perioperative evaluation of bleeding diathesis. Hematology Am Soc Hematol Educ Program 2006; 457-61.
3. Eikelboom JW, Hirsh J. Monitoring unfractionated heparin with the aPTT: Time for a fresh look. Thromb Haemost 2006; 96: 547-52
4. Bates SM, Johnston M, Hirsh J, Ginsberg JS. Use of a fixed activated partial thromboplastin time ratio to establish a therapeutic range for unfractionated heparin. Arch Intern Med 2001; 161: 385-91.
5. Hirsh J, Raschke R. Heparin and low-molecular-weight heparin. The seventh ACCP Conference on Antithrombotic and Thrombolytic Therapy. Chest 2004; 126: 188S-203S.
6. Olson JD, Arkin CF, Brandt JT, et al. College of American Pathologists Conference XXXI on laboratory monitoring of anticoagulant therapy: laboratory monitoring of unfractionated heparin therapy. Arch Pathol Lab Med. 1998; 122: 782-98
7. CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation Assays and Molecular Hemostasis Assays; Approved Guideline—Fifth Edition. CLSI document H21-A5 (ISBN 1-56238-657-3). CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2008.
8. Hattersley PG, Hayse D. The effect of increased contact activation time on the activated partial thromboplastin time. Am J Clin Pathol 1976; 66: 479-82.
9. Ambrus JL, Schimert G, Lajos TZ, et al. Effect of antifibrinolytic agents and estrogens on blood loss and blood coagulation factors during open heart surgery. J Med 1971; 2: 65-81.
10. Prasad RN, Koh SC, Viegas OA, Ratnam SS. Effects on hemostasis after two-year use of low dose combined oral contraceptives with gestodene or levonorgestrel. Clin Appl Thromb Hemost 1999; 51: 60-70.
11. Solomon GE, Hilgartner MW, Kutt H. Coagulation defects caused by diphenylhydantoin. Neurology 1972; 22: 1165-71.
12. Young DS, ed. Effects of Drugs on Clinical Laboratory Tests, 5th ed. Washington: AACC Press, 2000.
13. Greinacher A, Warkentin TE. The direct thrombin inhibitor hirudin. Thromb Haemost 2008; 99: 819-29.
14. Soloway HB, Cox SP, Donahoo JV. Sensitivity of the activated partial thromboplastin time to heparin: effect of anticoagulant used for sample collection. Am J Clin Pathol 1973; 59: 760-2.
15. Lippi G, Franchini M, Montagnana M et al. Quality and reliability of routine coagulation testing: can we trust that sample? Blood Coag Fibrinol 2006; 17: 513-9.
16. CLSI. One-Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline. CLSI document H47-A [ISBN 1-56238-301-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 1996.
17. CLSI. How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline - Second Edition. CLSI document C28-A2 [ISBN 1-56238-406-6]. CLSI, 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898, USA, 2000.
18. Mannucci PM, Canciani MT, Mari D, Meucci P. The varied sensitivity of partial thromboplastin and prothrombin time reagents in the demonstration of the lupus-like anticoagulant. Scand J Haematol. 1979; 22: 423-32.
19. Schleider MA, Nachman RL, Jaffe EA, Coleman M. A clinical study of the lupus anticoagulant. Blood 1976; 48: 499-509.

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