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- 1) Armstrong B, Zidar JP, Ohman EM: The use of intra-aortic balloon counterpulsation in acute myocardial infarction and high risk coronary angioplasty. *Journal Interventional Cardiology* 8:185,1995
- 2) Bauriedel G, Schwaiblmair M, Kreuzer E, et al: Percutaneous intraaortic counterpulsation as a therapeutic option in cardiogenic shock. *Dt Sch Med Wschr* 20:834, 1995
- 3) Cannon AD, Roubin GS, Iyer SS, Dean LS: Technical and clinical expertise in the use of coronary stents. *Cardio Intervention* 9:23, 1993
- 4) Dietl CA, Berkheimer MD, Woods EL, Gilbert CL, Pharr WF, Beniot CH: Efficacy and cost effectiveness of pre-operative IABP in patients with ejection fraction of 0.25 or less. *Annals of Thoracic Surgery* 62:401, 1996
- 5) Georgeson S, Tolbert-Coombs A, Eckman MH: Prophylactic use of the intra-aortic pump in high risk cancer patients undergoing non cardiac surgery. A decision analytic view. 6:92 665 *American Journal of Medicine* 1992
- 6) Goetz RH, Bregman D, Esrig B, Laniado S: Unidirectional intraaortic balloon pumping in cardiogenic shock and intractable left ventricular failure. *American Journal of Cardiology* 29:213, 1972
- 7) Gottlieb S: Intra-aortic balloon counterpulsation: An important component of community hospital care. *Cardiac Assists* 5:1, 1990
- 8) Grotz RL, Yeston NS: Intra-aortic balloon counterpulsation in high risk cardiac patients undergoing noncardiac surgery. *Surgery* 106:1.1, 1989
- 9) Gurbel PA, Anderson RD, MacCord CS, Scott H, Konjathy SF, Poulton J, Stafford JL, Godard J: Arterial diastolic pressure augmentation by intra-aortic balloon counterpulsation enhances the onset of coronary artery reperfusion by thrombolytic therapy. *Circulation* 89:1 361, 1994
- 10) Holmes DR, et al: Contemporary reperfusion therapy for cardiogenic shock: the Gusto-I Trial experience. *Journal of the American College of Cardiology* 26:668-674, 1995
- 11) Ishihara M, Sato H, Tateishi H, et al: Intra-aortic balloon pumping as adjunctive therapy to rescue coronary angioplasty after failed thrombolysis in anterior wall acute myocardial infarction. *American Journal Cardiology* 76:73, 1995
- 12) Iskandrian AS: Aortic Counterpulsation: Small and large balloons (editorial comment). *Catheterization and Cardiovascular Diagnosis* 36:88, 1995
- 13) Kantrowitz A, Tjonneland S, Freed PS, et al: Initial clinical experience with intra-aortic balloon pumping in cardiogenic shock. *JAMA* 203:135, 1968
- 14) Kaul U, Sahay S, Bahl VK, et al: Coronary angioplasty in high risk patients: Comparison of elective intra-aortic balloon pump and percutaneous cardiopulmonary bypass support - A randomized study. *J of Interventional Cardiology* 8:199, 1995
- 15) Kern MJ, Aquirre FV, Bach R, Donohue T, Siegel R, Segal J: Augmentation of coronary blood flow by intra-aortic balloon pumping in patients after coronary angioplasty. *Circulation* 81:2, 500, 1993
- 16) Kern MJ, Aquirre FV, Tatineni S, Penick D, Serota H, Donohue T, Walter K: Enhanced coronary blood flow velocity during intraaortic balloon counterpulsation in critically ill patients. *JACC* 21:2, 359, 1993
- 17) Kern MJ: Intra-aortic balloon counterpulsation. *Coronary Artery Disease* 2:6, 649, 1991
- 18) Mackenzie DJ, Wagner WH, Kulber DA, et al: Vascular complications of intra-aortic balloon pump. *Am J Surg* 164:517, 1992
- 19) Makhoul RG, Cole CW, McCann RL: Vascular complications of intra-aortic balloon pump: An analysis of 436 patients. *Am Surgeon* 59:564, 1993
- 20) McGee MG, Zillgit SL, Trono R et al: Retrospective analyses of the need for mechanical circulatory support after cardiopulmonary bypass: a 44 month study of 14, 168 patients. *American Journal of Cardiology* 46:135-142, 1980
- 21) Mercer D, Doris P, Salerno TA: Intra-aortic balloon counterpulsation in septic shock. *Can J Surg* 24:643, 1981
- 22) Ohman EM, George BS, White CJ: Use of aortic counterpulsation to improve sustained coronary artery patency during acute myocardial infarction. *Circulation* 90:2, 792, 1994
- 23) Ohman ME, Califf RM, George BS, et al: The use of intra-aortic balloon pumping as an adjunct to reperfusion therapy in acute myocardial infarction, *Am Heart J* (March):895-901, 1991
- 24) O' Murchu B, Foreman RD, Shaw RE, Brown DL, Petersen KL, Buchbinder M: Role of intra-aortic balloon pump counterpulsation in high risk coronary rotational atherectomy. *JACC* 26:5, 1270, 1995
- 25) Quaal SJ: *Comprehensive Intra-aortic Balloon Pumping*. ed 2, St Louis, MO: CV Mosby Co; 1993
- 26) Snow N, Lucas AE, Richardson JD: Intra-aortic balloon counterpulsation for cardiogenic shock from cardiac contusion. *The Journal at Trauma* 22:5, 426, 1982
- 27) Stomel RJ, Rasak M, Bates ER: Treatment strategies for acute myocardial infarction complicated by cardiogenic shock in a community hospital. *Chest* 105(4):997-1002, 1994
- 28) Underwood MJ, Firmin RK, Graham TR: Current concepts in the use of intra-aortic balloon counterpulsation. *BR J Hospital Med* 50:391, 1993
- 29) Waksman R, Weiss AT, Gotsman MS, Hasis Y: Intra-aortic balloon counterpulsation improves survival in cardiogenic shock complication acute myocardial infarction. *European Heart Journal* 14, 71, 1993
- 30) Reesink, K., Van der Nagel, T., Bovelander, J., Jansen, J., van der Veer, F, Schreuder, J.: Feasibility of a Fiber Optic System for Invasive Blood Pressure Measurements. *Catheterization and Cardiovascular Interventions*, 2002; 57: 272-6.
- 31) Wolthuis, R., Mitchell, G., Saaski, E., Hartl, J., Afromowitz, A.: Development of Medical Pressure and Temperature Sensors Employing Optical Spectrum Modulation. *IEEE Transactions on Biomedical Engineering*, 1991; 38(10): 974-81.
- 32) Schreuder, J., Donelli, A., Hanlon, P., Maisano, F., Castiglioni, A., Opizzi, M., Alfieri, O.: Closed Loop Intra-Aortic Balloon Counterpulsation in Patients with Marked Arrhythmia Using a Real Time Dicrotic Notch Prediction Algorithm. *The American Journal of Cardiology*, September 24, 2002; TCT Abstracts/Posters: 44H.
- 33) Donelli, A., Jansen, J., Hoeksel, B., Pedferri, P., Hanania, R., Bovelander, J., Maisano, F., Castiglioni, A., Alfieri, O., Schreuder, J.: Performance of a Real Time Dicrotic Notch Detection and Prediction Algorithm in Arrhythmic Human Aortic Pressure Signals. *Journal of Clinical Monitoring and Computing*, 2002; 17: 181-5

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Arrow International, Inc.