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- 3 Forbes, B.A., D.F. Sahm, and A.S Weissfeld. 1998. Bailey and Scott's diagnostic microbiology. 10th ed. Mosby, St. Louis, MO.
- 4 Wardford, A., M. Chernesky, and E. M. Peterson. 1999. Cumitech 19A, Laboratory diagnosis of Chlamydia trachomatis infections. ASM, Washington, DC.
- 5 Miller, J. M. 1999. A guide to specimen management in clinical microbiology, 2nd ed. ASM, Washington, DC.
- 6 Isenberg, H. D., 2004. Clinical microbiology procedures handbook, 2nd ed. ASM, Washington, DC.
- 7 Isenberg, H.D., 1998. Essential procedures for clinical microbiology. Chapter 14.12, Page 787. Packaging and shipping infectious substances.
- 8 Maass, M. and K. Dalhoff. 1995. Transport and storage conditions for cultural recovery of Chlamydia pneumoniae. J.Clin. Microbiol. 33:1793-1796.
- 9 Johnson, F. 2005. Transport of viral specimens. Clin. Microbiol. Rev 3:120-121.
- 10 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
- 11 National Committee for Clinical Laboratory Standards (NCCLS). 1994. Procedures for handling and transport of diagnostic specimens and etiologic agents. Approved Standard H5-A3, NCCLS, Wayne, PA..
- 12 National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality control of microbiological transport systems. Approved Standard M40-A, NCCLS, Wayne, PA.
- 13 National Committee for Clinical Laboratory Standards (NCCLS). 2004. Viral culture. Proposed Standard M41, NCCLS, Wayne, PA.

BD, BD Logo, BBL, and CultureSwab are trademarks of Becton, Dickinson and Company.©2002 BD.
ATCC is a trademark of the American Type Culture Collection.

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