





## Precautions when handling insulin syringes

(Case 4) A medical staff member usually administered mixed insulin measured up to the maximum dose of a 30-unit insulin syringe. However, on that day, the staff member erroneously took a 50-unit insulin syringe and mixed the maximum dose as usual.

## Make sure to confirm the size (unit) of insulin syringes.



Insulin syringes are delivered in several sizes, each with different maximum doses of the unit. Review the unit of insulin syringes to be used to prevent a mix-up. For example, by unifying the unit of insulin syringes used.

The Ministry of Health, Labour and Welfare (MHLW) issued a notification related to PMDA Medical Safety Information No. 23:

●PSEHB/PSD Notification No. 0519-1 dated on May 19, 2020 Revision of Precautions

## About this information

- PMDA Medical Safety Information is issued by the Pharmaceuticals and Medical Devices Agency for the purpose of providing healthcare providers with clearer information from the perspective of promoting the safe use of pharmaceuticals and medical devices. The information presented here has been compiled, with the assistance of expert advice, from cases collected as Medical Accident Information Reports by the Japan Council for Quality Health Care, and collected as Adverse Drug Reaction and Malfunction Reports in accordance with the Law on Securing Quality, Efficacy and Safety of Pharmaceuticals and Medical Devices.
- \* We have tried to ensure the accuracy of this information at the time of its compilation but do not guarantee its accuracy in the future.
- \* This information is not intended to impose constraints on the discretion of healthcare professionals or to impose obligations and responsibility on them, but is provided as a support to promote the safe use of pharmaceuticals and medical devices by healthcare professionals.

Access to the most up-to-date safety information is provided via the PMDA Medinavi service.



Published by the Pharmaceuticals and Medical Devices Agency

Contact: Medical Safety Information Group

4/4

TEL 03-3506-9486 FAX 03-3506-9543 https://www.pmda.go.jp/