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○ NCI(National Cancer Institute) Fact sheet: Cryosurgery in Cancer Treatment: Questions and Answers

「Cryosurgery may offer an option for treating cancers that are considered inoperable or that do not respond to standard treatments. Furthermore, it can be used for patients who are not good candidates for conventional surgery because of their age or other medical conditions.」

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○ NCCN(National Comprehensive Cancer Network) Clinical Practice Guidelines in Oncology

「a small set of elderly or infirm patients with small tumors may be offered surveillance alone or energy ablative techniques, such as radiofrequency ablation or cryoablation.」

○EAU(European association of Urology) EAU Guideline

「Currently, patients not suitable for open or laparoscopic surgery due to poor performance status with smaller peripheral tumors should be considered for the above-mentioned [Image-guided percutaneous and minimally invasive techniques, e.g. percutaneous radiofrequency (RF) ablation, cryoablation, microwave ablation, laser ablation and high-intensity focused ultrasound ablation (HIFU)] techniques for RCC treatment (grade B recommendation [Based on well-conducted clinical studies, but without randomized clinical trials]).」

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- 1) Seifert JK, Morris DL: World survey on the complications of hepatic and prostate cryotherapy, *World Journal of Surgery*, 23(2)109-113, 1999
- 2) Sarantou T, Bilchik A, Ramming KP: Complications of hepatic cryosurgery, *Seminars in Surgical Oncology*, 14(2)156-162, 1998
- 3) Carvalhal EF, Gill IS, Meraney AM, et al: Laparoscopic renal cryoablation: impact on renal function and blood pressure, *Urology*, 58(3)357-361, 2001
- 4) Seifert JK, Junginger T, Morris DL: A collective review of the world literature on hepatic cryotherapy; *Journal of the Royal College of Surgeons of Edinburgh*; 43 (3) 141-154; 1998
- 5) Shingleton WB, Sewell PE Jr: Percutaneous renal tumor cryoablation with magnetic resonance imaging guidance; *The Journal of Urology*; 165 (3) 773-776; 2001
- 6) Cowan BD, Sewell PE, Howard JC, et al: Interventional magnetic resonance imaging cryotherapy of uterine fibroid tumors: Preliminary observation; *American Journal of Obstetrics and Gynecology*; 186 (6) 1183-118; 2002

- 7) Warlick CA, Lima GC, Allaf ME, et al; Clinical sequelae of radiographic iceball involvement of collecting system during computed tomography-guided percutaneous renal tumor cryoablation; *Urology*; 67 (5) 918-922; 2006
- 8) 岸本幸一、波多野孝史、三木健太、他; 腎腫瘍に対するMRガイド下凍結療法の腎機能に関する検討; 第33回日本低温医学会総会抄録; p95; 2006
- 9) Rukstalis DB, Khorsandi M, Garcia FU, et al; Clinical experience with open renal cryoablation; *Urology*; 57 (1) 34-39; 2001
- 10) Gill IS, Novick AC, Meraney AM, et al; Laparoscopic renal cryoablation in 32 patients; *Urology*; 56 (5) 748-753; 2000
- 11) Lee DI, McGinnis DE, Feld R, et al; Retroperitoneal laparoscopic cryoablation of small renal tumors: intermediate results; *Urology*; 61 (1) 83-88; 2003
- 12) Shingleton WB, Louisville KY, Sewell PE, et al; Percutaneous renal tumor cryoablation: results in the first 90 patients; *The journal of urology*; 171(Suppl)463, A1751; 2004.
- 13) Gill IS, Remer EM, Hasan WA, et al; Renal cryoablation: outcome at 3 years; *The journal of Urology*; 173(6)1903-1907; 2005
- 14) Cestari A, Guazzoni G, Dell'acqua V, et al; Laparoscopic cryoablation of solid renal masses: intermediate term followup; *The journal of urology*; 172(4 Pt 1)1267-1270; 2004
- 15) Silverman SG, Tuncali K, van Sonnenberg E, et al; Renal tumors: MR Imaging-guided percutaneous cryotherapy – Initial experience in 23 patients; *Radiology*; 236(2)716-724; 2005
- 16) Khorsandi M, Foy RC, Chong W, et al; Preliminary experience with cryoablation of renal lesions smaller than 4 centimeters; *Journal of the american osteopathic association*; 102(5)277-281; 2002
- 17) Moon TD, Lee FT, Hedican SP, et al; Laparoscopic cryoablation under sonographic guidance for the treatment of small renal tumors; *Journal of Endourology*; 18(5)436-440; 2004
- 18) Nadler RB, Kim SC, Rubenstein JN, et al; Laparoscopic renal cryosurgery: the northwestern experience; *The journal of Urology*; 170(4 Pt 1)1121-1125; 2003
- 19) Long JP, Faller GT; Percutaneous Cryoablation of the Kidney in a Porcine Model; *Cryobiology*; 38 (1) 89-93; 1999
- 20) Shingleton WB, Farabaugh P, Hughson M et al; Effects of cryoablation on short-term development of urinary fistulas in the porcine kidney.; *Journal of Endourology*; 17 (1) 37-40; 2003
- 21) Remer EM, Sung GT, Meraney AM et al; Effect of intentional cryoinjury to the

- renal collecting system; *The Journal of Urology*; 163 (suppl 4); 2000
- 22) Ladd AP, Rescorla FJ, Baust JG et al; Cryosurgical effects on growing vessels; *The American Surgeon*; 65 (7) 677-682; 1999
 - 23) Sewell P, Shingleton W, Five-year treatment success and survival of patients treated with percutaneous IMRI guided and monitored renal cell carcinoma Cryoablation. *BJU Int.* 94:106. Abstract MP-1609, 2004
 - 24) Hegarty, N J, Gill, I S, Desai, M M, Remer, E M, O'Malley, C M, Kaouk, J H: Probe-ablative nephron-sparing surgery: cryoablation versus radiofrequency ablation.. *Urology* 68:7-13,2006
 - 25) Hafron J, Kaouk JH. Cryosurgical ablation of renal cell carcinoma.*Cancer Control.* 2007 Jul;14(3):211-7.
 - 26) Atwell TD, Farrell MA, Leibovich BC, Callstrom MR, Chow GK, Blute ML, Charboneau JW: Percutaneous renal cryoablation: experience treating 115 tumors.. *J Urol* 179:2136-2140,2008
 - 27) Davol, P E, Fulmer, B R, Rukstalis, D B: Long-term results of cryoablation for renal cancer and complex renal masses.. *Urology* 68:2-6,2006
 - 28) Finley, D S, Beck, S: Percutaneous and laparoscopic cryoablation of small renal masses.. *J Urol* 180:492-498,2008
 - 29) Bandi, G, Wen, C C, Hedican, S P, Moon, T D, Lee, F T Jr, Nakada, S Y: Cryoablation of small renal masses: assessment of the outcome at one institution.. *BJU Int* 100:798-801,2007
 - 30) Lawatsch, E J, Langenstroer, P, Byrd, G F, See, W A, Quiroz, F A, Begun, F P: Intermediate results of laparoscopic cryoablation in 59 patients at the Medical College of Wisconsin.. *J Urol* 175:1225-1229,2006
 - 31) Johnson DB, Solomon SB, Su LM, Matsumoto ED, Kavoussi LR, Nakada SY, Moon TD, Shingleton WB, Cadeddu JA. Defining the complications of cryoablation and radio frequency ablation of small renal tumors: a multi-institutional review.. *J Urol* 172:874-877, 2004
 - 32) Bandi G, Hedican SP, Nakada SY: Current practice patterns in the use of ablation technology for the management of small renal masses at academic centers in the United States.. *Urology* 71:113-117,2008
 - 33) Pantuck AJ, Zisman A, Cohen J, Belldegrun A. Cryosurgical ablation of renal tumors using 1.5-millimeter, ultrathin cryoprobes. *Urology.* 2002 Jan;59(1):130-3.
 - 34) Weld KJ, Figenshau RS, Venkatesh R, Bhayani SB, Ames CD, Clayman RV, Landman J. Laparoscopic cryoablation for small renal masses: three-year follow-up. *Urology.* 2007 Mar;69(3):448-51.

- 35) Wyler SF, Sulser T, Ruzsat R, Weltzien B, Forster TH, Provenzano M, Gasser TC, Bachmann A. Intermediate-term results of retroperitoneoscopy-assisted cryotherapy for small renal tumours using multiple ultrathin cryoprobes. *Eur Urol.* 2007 Apr;51(4):971-9. Epub 2006 Oct 27.
- 36) Bachmann A, Sulser T, Jayet C, Wyler S, Ruzsat R, Reich O, Gasser TC, Siebels M, Stief CG, Casella R. Retroperitoneoscopy-assisted cryoablation of renal tumors using multiple 1.5 mm ultrathin cryoprobes: a preliminary report. *Eur Urol.* 2005 Apr;47(4):474-9. Epub 2005 Jan 13.
- 37) Morrison PR, Silverman SG, Tuncali K, Tatli S. MRI-guided cryotherapy. *J Magn Reson Imaging.* 2008 Feb;27(2):410-20. Review.