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- 1 Guidance for the Non-Clinical and Clinical Investigation of Devices Used for the Treatment of Benign Prostatic Hyperplasia (BPH). FDA Guidance issued on August 17, 2010.
- 2 National Institute of Diabetes and Digestive and Kidney Diseases
<https://www.niddk.nih.gov/health-information/urologic-diseases/prostate-problems/prostate-enlargement-benign-prostatic-hyperplasia>
- 3 Roehrborn CG. Benign prostatic hyperplasia: an overview. Rev Urol. 2005;7 Suppl 9(Suppl 9):S3-S14.
- 4 Yukio Homma et al. Epidemiologic Survey of lower Urinary Tract Symptoms in Asia and Australia Using the International Prostate Symptom Score: Int J Urol 1997;4:40-46
- 5 Clinical Summary Cantrill Rezum Post-Procedure Care URO-724504-AA
- 6 James U, Melissa M et al. Cost-effectiveness analysis of six therapies for the treatment of lower urinary tract symptoms due to benign prostatic hyperplasia. ClinicoEconomics and Outcomes Research 2018;10 29-43
- 7 J. Kellogg Parsons. Benign Prostatic Hyperplasia and Male Lower Urinary Tract Symptoms: Epidemiology and Risk Factors. Curr Bladder Dysfunct Rep (2010) 5: 212-218
- 8 Shaun W Lee, Esther M Chan et al. The global burden of lower urinary tract symptoms suggestive of benign prostatic hyperplasia: A systematic review and meta-analysis. Scientific Reports 7: 7984
- 9 厚生労働省：平成 29 年患者調査（傷病分類編）
- 10 男性下部尿路症状・前立腺肥大症診療ガイドライン
- 11 Gacci M, Corona G et al. Metabolic syndrome and benign prostatic enlargement: a systematic review and meta-analysis. BJU Int 2015; 115: 24–31
- 12 Masumori N, Tsukamoto T et al. Japanese men have smaller prostate volumes but comparable urinary flow rates relative to American men: results of community-based studies in 2 countries. J Urol 1996; 155: 1324–1327

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著者	重みづけ平均が計算されている評価項目	重みづけ平均の計算に必要な情報が記載されている評価項目
Marra ら	Retrograde ejaculation rate	-
Cornu ら	-	B-TURP vs M-TURP; (A) interventional duration, min); (B) transfusion rate; (C) hemoglobin loss; (D) sodium decrease; (E) urinary tract infections; (F) immediate acute urinary retention; (G) clot retention; (H) catheterization time; (I) recatheterization; (J) immediate reoperation rate; (K) transurethral resection syndrome; (L) length of stay. (M - Q) Efficacy at 12 mo: (M) International Prostate Symptom Score (IPSS) at 12 mo; (N) maximum flow rate (Qmax) at 12 mo; (O) quality of life at 12 mo; (P) prostate volume at 12 mo; (Q) postvoid residual (PVR) at 12 mo. (R-T) Complications at 12 mo: (R) strictures at 12 mo; (S) incontinence rate at 12 mo; (T) reoperation at 12 mo. (U-W) Long-term efficacy: (U) long-term IPSS; (V) long-term Qmax; (W) long-term PVR. (X-Z) Long-term complications: (X) strictures; (Y) bladder neck contracture; (Z) reoperation rate. HoLEP vs M-TURP (short-term complications); (A) Intervention duration; (B) sodium decrease; (C) hemoglobin decrease; (D) acute urinary retention; (E) clot retention; (F) recatheterization; (G) reoperation; (H) catheter duration; (I) length of stay; (J) storage symptoms; (K) strictures after 1 yr; (L) urinary tract infections. HoLEP vs M-TURP (Efficacy); (A) International Prostate Symptom Score (IPSS) difference after 1-yr follow-up; (B) maximum flow rate (Qmax) difference after 1-yr follow-up; (C) postvoid residual difference after 1-yr follow-up; (D) quality of life difference after 1-yr follow-up; (E) IPSS difference after long-term follow-up; (F) Qmax difference after long-term follow-up. PVP vs M-TURP; (A) Intervention duration; (B) transfusion rate; (C) urinary tract infection; (D) acute urinary retention; (E) catheter duration; (F) length of stay; (G) bladder neck contracture (long term); (H) stricture (long-term rate); (I) reoperation for bladder outlet obstruction relief (long-term rate); (J) International Prostate Symptom Score at 12 mo; (K) maximum flow rate at 12 mo; (L) postvoid residual at 12 mo.
Laurenco ら	-	B-TURP vs M-TURP; Symptom score at 12 months, Peak urine flow rate (ml/s) at 12 months HoLEP vs M-TURP; Symptom score at 12 months, Peak urine flow rate (ml/s) at 12 months PVP vs M-TURP Symptom score at 12 months, Peak urine flow rate (ml/s) at 12 months

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4. Marra G, Sturch P, Oderda M, Tabatabaei S, Muir G, Gontero P. Systematic review of lower urinary tract symptoms/benign prostatic hyperplasia surgical treatments on men' s ejaculatory function: time for a bespoke approach? *Int J Urol* 2016; 23: 22- 35
5. Lourenco, T., et al. Alternative approaches to endoscopic ablation for benign enlargement of the prostate: systematic review of randomised controlled trials. *Bmj*, 2008. 337: a449.
6. Cornu, J.N., et al. A Systematic Review and Meta-analysis of Functional Outcomes and Complications Following Transurethral Procedures for Lower Urinary Tract Symptoms Resulting from Benign Prostatic Obstruction: An Update. *Eur Urol*, 2015. 67: 1066.

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