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- 1 National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- 2 Garner, J.S. 1996. Hospital Infection Control Practices Advisory

Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. *Infect. Control Hospital Epidemiol.* 17: 53-80.

- 3 U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- 4 Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). *Official Journal L262*, 17/10/2000, p. 0021-0045.
- 5 Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover. 1999. *Manual of clinical microbiology*. 7th ed. ASM, Washington, D.C.
- 6 Gleaves, C.A., R.L. Hodinka, S.L.G. Johnston, and E.M. Swierkosz. 1994. Cumitech 15A. Laboratory diagnosis of viral infections. ASM, Washington, DC.
- 7 Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 1998. *Bailey and Scott's diagnostic microbiology*. 10th ed. Mosby, St. Louis, MO.
- 8 Wardford, A., M. Chernesky, and E. M. Peterson. 1999. Cumitech 19A, Laboratory diagnosis of *Chlamydia trachomatis* infections. ASM, Washington, DC.
- 9 Miller, J. M. 1999. A guide to specimen management in clinical microbiology, 2nd ed. ASM, Washington, DC.
- 10 Isenberg, H. D., 2004. *Clinical microbiology procedures handbook*, 2nd ed. ASM, Washington, DC.
- 11 Isenberg, H.D., 1998. *Essential procedures for clinical microbiology*. Chapter 14.12, Page 787. Packaging and shipping infectious substances.
- 12 Maass, M. and K. Dalhoff. 1995. Transport and storage conditions for cultural recovery of *Chlamydia pneumoniae*. *J. Clin. Microbiol.* 33:1793-1796.
- 13 Johnson, F. 2005. Transport of viral specimens. *Clin. Microbiol. Rev* 3:120-121.
- 14 42CFR72. Code of Federal Regulations, Title 42, Volume 1, Part 72. Interstate Shipment of Etiologic Agents.
- 15 National Committee for Clinical Laboratory Standards (NCCLS). 1994. Procedures for handling and transport of diagnostic specimens and etiologic agents. Approved Standard H5-A3, NCCLS, Wayne, PA..
- 16 National Committee for Clinical Laboratory Standards (NCCLS). 2003. Quality control of microbiological transport systems. Approved Standard M40-A, NCCLS, Wayne, PA.
- 17 National Committee for Clinical Laboratory Standards (NCCLS). 2004. Viral culture. Proposed Standard M41, NCCLS, Wayne, PA.

BD, BD Logo, BBL, and CultureSwab are trademarks of Becton, Dickinson and Company. ©2002 BD.

ATCC is a trademark of the American Type Culture Collection.

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(COPAN Italia, S.p.A)