


# Clinical practice differences in Japan and USA when using open surgical grafts

HBD East 2023

December 14<sup>th</sup>

Ariake Central Tower Hall & Conference

# Synthetically Sealed Graft (SSG) – Key Benefits



**First new Open Surgical Graft Technology in Decades**

# Clinical Practice Similarities

---

- Same Patient Population
- Same disease states / conditions
- Same technique / procedure
- Same OR support staff
- Same sutures / suturing technique

**Therefore.... No additional training required**

# Clinical Practice Differences with SSG vs. Traditional Open Surgical Grafts



---

- Less Blood Loss
- Fewer Blood Products per Procedure
- Potential for Lower Risk of Infections
- Shorter on Pump and Cross Clamp times
- Shorter Procedure Times Free up Valuable OR space
- Less Traumatic Revisions if Required in the Future





# Suture Hole Leak Test Data

# Experimental Design

- A suture was sewn around each sample to mimic an end-to-end anastomosis
  - Average of 108 holes/sample
- 3 different sutures/needle sizes were used
  - 3-0, 4-0, & 5-0
- The grafts were then pressurized to 120 mm Hg for 1 minute with blood analog
  - 60/40 water/glycerol
- The leakage was captured and measured




# Results



# Conclusion

- If a typical anastomosis requires **20 suture holes** with a **4-0** suture and takes **7 minutes** to reach hemostasis...





# Impact of Reduced Blood Loss

---

- Patient benefits
  - Typical risks of bleeding during surgery can include complications associated with a higher risk of infection, ischemic events attributable to hypo-perfusion (e.g., myocardial infarction, acute kidney injury), in-hospital mortality, and transfusion-related adverse events. Bleeding can lead to cardiac tamponade, shock, dilutional coagulopathy and impaired hemostasis.
- Hospital benefits
  - Shorter time to hemostasis may lead to more cases performed per year
  - Shorter time to patient discharge, i.e. less DRG used on follow up care
  - Reduced cost / usage of blood products and surgical adhesives
- Economics of blood products and surgical adhesives
  - See U of FLA analysis

# Cost of Blood Products and Surgical Adhesives

Less time to hemostasis = potentially ↓ OR costs, ↑ cases per OR

**Average Cost: Chronic & Acute vs. Subacute/Chronic vs. Acute/Hyperacute/Acute on Chronic**

|  | Cyroprecipitate | Platelets  | Red Blood Cells | Fresh Frozen Plasma | Crossmatch | Other   | Blood Products (B-G) | Blood Products Only (B-E) | Surgical Adhesives | Clotting Factors | Grand Total: Blood Products + Surgical Adhesives + Clotting Factors |
|--|-----------------|------------|-----------------|---------------------|------------|---------|----------------------|---------------------------|--------------------|------------------|---|
| <b>Chronic &amp; Acute</b>               | \$2,269.83      | \$1,477.97 | \$1,251.59      | \$57.74             | \$372.46   | \$52.84 | \$5,274.25           | \$5,057.12                | \$1,048.43         |                  | <b>\$6,530.47</b>   |
| <b>Subacute/Chronic</b>                  | \$2,353.89      | \$1,532.71 | \$1,297.95      | \$59.87             | \$384.90   | \$54.81 | \$5,357.39           | \$5,244.42                | \$1,069.69         |                  | <b>\$6,687.55</b>   |
| <b>Acute/Hyperacute/Acute on Chronic</b> | \$2,278.96      | \$1,488.60 | \$1,253.49      | \$55.86             | \$372.23   | \$53.22 | \$5,292.70           | \$5,076.92                | \$1,050.58         |                  | <b>\$6,552.56</b>   |

**Average Cost: TEVAR Subacute/Chronic vs. TEVAR Acute/Hyperacute/Acute on Chronic**

|  | Cyroprecipitate | Platelets  | Red Blood Cells | Fresh Frozen Plasma | Crossmatch | Other   | Blood Products (B-G) | Blood Products Only (B-E) | Surgical Adhesives | Clotting Factors | Grand Total: Blood Products + Surgical Adhesives + Clotting Factors |
|--|-----------------|------------|-----------------|---------------------|------------|---------|----------------------|---------------------------|--------------------|------------------|---|
| <b>Subacute/Chronic</b>                  | \$2,648.38      | \$1,570.59 | \$1,177.63      | \$25.72             | \$345.77   | \$47.55 | \$5,412.51           | \$5,422.33                | \$1,062.43         |                  | <b>\$6,795.30</b>   |
| <b>Acute/Hyperacute/Acute on Chronic</b> | \$2,301.29      | \$1,501.50 | \$1,262.78      | \$54.93             | \$368.80   | \$54.01 | \$5,330.58           | \$5,120.50                | \$1,047.67         |                  | <b>\$6,590.59</b>   |

**Average Cost: TEVAR vs. Non TEVAR for Chronic and Acute Dissections**

|                                    | Cyroprecipitate | Platelets  | Red Blood Cells | Fresh Frozen Plasma | Crossmatch | Other   | Blood Products (B-G) | Blood Products Only (B-E) | Surgical Adhesives | Clotting Factors | Grand Total: Blood Products + Surgical Adhesives + Clotting Factors |
|------------------------------------|-----------------|------------|-----------------|---------------------|------------|---------|----------------------|---------------------------|--------------------|------------------|---|
| <b>Chronic and Acute w TEVAR</b>   | \$2,301.29      | \$1,501.50 | \$1,262.78      | \$54.93             | \$368.80   | \$54.01 | \$5,330.58           | \$5,120.50                | \$1,047.67         |                  | <b>\$6,590.59</b>   |
| <b>Chronic and Acute w/o TEVAR</b> | \$2,371.55      | \$1,530.65 | \$1,285.92      | \$61.23             | \$381.11   | \$55.57 | \$5,246.37           | \$5,249.35                | \$1,062.15         |                  | <b>\$6,552.49</b>   |

\* NOTE – Clotting Factors are a significant cost to any procedure when used

# Synthetically Sealed Graft – Infection Control

---


## Infection Control – Initial Testing

### Evaluation of the Anti-biofilm Property of Diaxamed Synthetically Sealed Grafts

The purpose of the study is to evaluate the anti-biofilm formation properties of Diaxamed Grafts and their ability to prevent bacterial migration from the outer surface to the inner surface of the grafts

# Anti-biofilm Property of SSG

1. Diaxamed and the competitor's grafts of a diameter of 6 mm were evaluated in this study
2. Individual grafts submerged in sterile 50 mL conical tubes containing 10 mL of growth media containing 1,000 CFU/mL overnight culture of *Staph. aureus*
3. At 24 hours post-incubation at 37°C, the bacterial culture in the sample reservoir was discarded, and the grafts were rinsed 3 times with 10 mL of phosphate-buffered saline (PBS) and gently mixed to remove any unbound cells
4. Bacteria biofilm on the outer surface of the washed grafts was recovered by adding 10 mL of PBS into the conical tubes and sonicated in an ultrasound water bath for 15 minutes, then enumerated
5. The experiment was conducted in triplicates independently, and the number of bacteria recorded




Diaxamed  
SSG

Competitor

# Bacteria Migration Prevention Property of SSG

1. Bacteria, if any, that managed to migrate from the outer surface of the grafts to the inner surface of the grafts were enumerated by adding 4 mL of PBS into individual grafts and sonicated in an ultrasound water bath for 15 minutes, then enumerated
2. The experiment was conducted in triplicates independently, and **the results indicated that no bacteria were able to migrate through the Diaxamed grafts**, i.e., no bacteria were recovered from the inside of the Diaxamed grafts




Diaxamed  
SSG

Competitor

# Anti-biofilm Testing Results


## Anti-biofilm Property of Diaxamed Synthetically Sealed Grafts



The bacteria recovered from the outer surface of the competitor's grafts are  $3.32 \times 10^6$  CFU/Cm<sup>2</sup>. The bacteria recovered from the outer surface of the Diaxamed's grafts are  $1.61 \times 10^6$  CFU/Cm<sup>2</sup>, which is **SSG has 52% less bacteria on the outer surface compared to Gelweave**.

# Bacterial Migration Prevention Results

## Bacterial Migration Prevention Property of Diaxamed Synthetically Sealed Grafts



The bacteria recovered from the inner surface of the competitor's grafts are  $6.87 \times 10^2$  CFU/Cm<sup>2</sup>.

**No bacteria were recovered from the inner surface of the Diaxamed's grafts, which indicates that the Diaxamed Synthetically Sealed Grafts were able to prevent the migration of bacteria from the outer surface to the inner surface**

# Interim GLP Pathology Report Summary

---




## **A Preclinical Study for the Chronic Evaluation of a New Thoracic Graft in the Sheep Model 180 day results**

Confidential Information – All information provided under NDA and not to be disclosed to any third party



# Synthetically Sealed Graft – GLP Study Design

- Sheep aortic model
- Sheep 1-3 years / 75kg-95kg
- End to end anastomosis
- 12x test (SSG, woven conduit)
- 5x control (Gelweave, woven conduit)
- GLP Site – University of Minnesota



# Identification Information for all 180-day samples

| University of Minnesota ID | CVP Accession Number | Target Tissue | Type of Graft Implanted | Time Point/ Control (C) or Test (T) |
|----------------------------|----------------------|---------------|-------------------------|-------------------------------------|
| <b>DTG-1</b>               | DXM23-009            | Aortic Graft  | Gelweave 18mm           | 180 Day (Control)                   |
| <b>DTG-2</b>               | DXM23-010            | Aortic Graft  | Gelweave 18mm           | 180 Day (Control)                   |
| <b>DTG-3</b>               | DXM23-011            | Aortic Graft  | SSG 16mm                | 180 Day (Test)                      |
| <b>DTG-4</b>               | DXM23-012            | Aortic Graft  | SSG 20mm                | 180 Day (Test)                      |
| <b>DTG-6</b>               | DXM23-013            | Aortic Graft  | Gelweave 18mm           | 180 Day (Control)                   |
| <b>DTG-7</b>               | DXM23-014            | Aortic Graft  | SSG 20mm                | 180 Day (Test)                      |
| <b>DTG-8</b>               | DXM23-015            | Aortic Graft  | SSG 20mm                | 180 Day (Test)                      |
| <b>DTG-9</b>               | DXM23-016            | Aortic Graft  | SSG 20mm                | 180 Day (Test)                      |
| <b>DTG-10</b>              | DXM23-008            | Aortic Graft  | SSG 20mm                | *150 Day (Test)                     |
| <b>DTG-11</b>              | DXM23-017            | Aortic Graft  | SSG 20mm                | 180 Day (Test)                      |
| <b>DTG-13</b>              | DXM23-018            | Aortic Graft  | SSG 20mm                | 180 Day (Test)                      |


\* Cause of early death is not device related, final determination of cause to be made however suspected to be confirmed adrenal disorder.

## In Situ Graft Images - At Implantation

---




SSG- No graft Blushing, entirely blood tight



Gelweave - Graft Blushing

# Summary of On Pump and Cross Clamp Time all implants

## Controls Vs Test On Pump and Cross Clamp Durations



| Category           | # of implants | Average Cross clamp time (minutes) | Average On Pump time (minutes) |
|--------------------|---------------|------------------------------------|--------------------------------|
| Control (Gelweave) | 5             | 46.8                               | 78.32                          |
| Test (SSG)         | 12            | 40.17                              | 63.92                          |
| % Time Difference  | -             | 16.5% Increase                     | 22.5% Increase                 |

## ***Along with Clinical and Economic benefits, Future benefits for Doctors in both Japan and USA***

---

- “Nice , not as stiff as expected , easy to cut and size, like the orientation line along the graft “
- “No dulling of the needle. Usually by second anastomosis, the drag on the needle is much greater. This is not the case with the SSG.”
- “When suturing and you release the needle to grab it from the other side, it does not move or fall out. Gives you the ability to have a consistent and effortless cadence.”
- “Saturability and handling very nice”
- “Nice, easy to cut, sew and handle in the surgical field”
- “Graft appeared hemostatic instantly when clamps were removed”

Implantation Notes on Diaxamed SSG During GLP Study

# Thank You!