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1. Clinical and Laboratory Standards Institute (formerly NCCLS). *Interference Testing in Clinical Chemistry: Approved Guideline*. NCCLS document EP7-A. Villanova, PA: NCCLS, 2002.
2. Lott JA, Abbott LB. Creatine Kinase Isoenzymes. *Clin Lab Med* 1986 Sept;6(3):547-76.
3. Clinical and Laboratory Standards Institute (formerly NCCLS). *How to Define and Determine Reference Intervals in the Clinical Laboratory; Approved Guideline – Second Edition*. NCCLS Document C28-A2. Wayne, PA: NCCLS, 2000.
4. Galen RS and Gambino SR. Isoenzymes of CPK and LDH in myocardial infarction and certain other diseases. *Pathobiol Annu* 1975; 5:283-315.
5. Gibler WB, Young GP, Hedges JR, et al. Acute myocardial infarction in chest pain patients with nondiagnostic ECGs: serial CK-MB sampling in the emergency department. *Ann Emerg Med* 1992;21(5):504-12.
6. Hedges JR, Young GP, Henkel GF, et al. Serial ECGs are less accurate than serial CK-MB results for emergency department diagnosis of myocardial infarction. *Ann Emerg Med* 1992;21(12):1445-50.
7. Apple FS. Acute myocardial infarction and coronary reperfusion. Serum cardiac markers for the 1990s. *Am J Clin Pathol* 1992;97(2):217-26.
8. Garabedian HD, Gold HK, Yasuda T, et al. Detection of coronary artery reperfusion with creatine kinase-MB determinations during thrombolytic therapy: correlation with acute angiography. *J Am Coll Cardiol* 1988;11(4):729-34.
9. Lewis BS, Ganz W, Laramie P, et al. Usefulness of a rapid initial increase in plasma creatine kinase activity as a marker of reperfusion during thrombolytic therapy for acute myocardial infarction. *Am J Cardiol* 1988;62:20-4.
10. Jockers-Wretou E, Pfeleiderer G. Quantitation of creatine kinase isoenzymes in human tissues and sera by an immunological method. *Clin Chim Acta* 1975;58:223-32.
11. Ogunro EA, Hearse DJ, Shillingford JP. Creatine kinase isoenzymes: their separation and quantitation. *Cardiovasc Res* 1977;11:94-102.
12. Tsung SH. Creatine kinase isoenzyme patterns in human tissue obtained at surgery. *Clin Chem* 1976;22:173-5.
13. Tsung JS, Tsung SS. Creatine kinase isoenzymes in extracts of various human skeletal muscles. *Clin Chem* 1986;32(8):1568-70.
14. Apple FS, Rogers MA, Sherman WM, et al. Profile of creatine kinase isoenzymes in skeletal muscles of marathon runners. *Clin Chem* 1984;30(3):413-6.
15. Ohman EM, Teo KK, Johnson AH, et al. Abnormal cardiac enzyme responses after strenuous exercise: alternative diagnostic aids. *Br Med J* 1982;285(6354):1523-6.

16. Apple FS, Rogers MA, Sherman WM, et al. Comparison of serum creatine kinase and creatine kinase MB activities post marathon race versus post myocardial infarction. *Clin Chim Acta* 1984;138:111-8.
17. Jaffe AS, Garfinkel BT, Ritter CS, et al. Plasma MB creatine kinase after vigorous exercise in professional athletes. *Am J Cardiol* 1984;53:856-8.
18. Wu AHB, Wang XM, Gornet TG, et al. Creatine kinase MB isoforms in patients with skeletal muscle injury: ramifications for early detection of acute myocardial infarction. *Clin Chem* 1992;38(12):2396-400.
19. Larca LJ, Coppola JT, Honig S. Creatine kinase MB isoenzyme in dermatomyositis: a noncardiac source. *Ann Int Med* 1981;94:341-3.
20. Keshgegian AA (editor). Serum creatine kinase MB isoenzyme in chronic muscle disease. *Clin Chem* 1984;30(4):575-8.
21. Yasminieh WG, Ibrahim GA, Abbasnezhad M, et al. Isoenzyme distribution of creatine kinase and lactate dehydrogenase in serum and skeletal muscle in duchenne muscular dystrophy, collagen disease, and other muscular disorders. *Clin Chem* 1978;24(11):1985-9.
22. Jaffe AS, Ritter C, Meltzer V, et al. Unmasking artifactual increases in creatine kinase isoenzymes in patients with renal failure. *J Lab Clin Med* 1984;104(2):193-202.
23. Tsung SH. Several conditions causing elevation of serum CK-MB and CK-BB. *Am J Clin Pathol* 1981;75(5):711-5.
24. Kimler SC, Sandhu RS. Circulating CK-MB and CK-BB isoenzymes after prostate resection. *Clin Chem* 1980;26(1):55-9.
25. Potkin RT, Werner JA, Trobaugh GB, et al. Evaluation of noninvasive tests of cardiac damage in suspected cardiac contusion. *Circulation* 1982;66(3):627-31.
26. El Allaf M, Chapelle JP, EL Allaf D, et al. Differentiating muscle damage from myocardial injury by means of the serum creatine kinase (CK) isoenzyme MB mass measurement/total CK activity ratio. *Clin Chem* 1986;32(2):291-5.
27. Schroff RW, Foon KA, Beatty SM, et al. Human anti-murine immunoglobulin responses in patients receiving monoclonal antibody therapy. *Cancer Res* 1985;45:879-85.
28. Boscato LM and Stuart MC. Heterophilic antibodies; a problem for all immunoassays. *Clin Chem* 1988;34(1):27-33.
29. Clinical and Laboratory Standards Institute (formerly NCCLS). *Evaluation of Precision Performance of Clinical Chemistry Devices - Approved Guideline*. NCCLS document EP5-A. Wayne, PA: NCCLS, 1999.
30. Panteghini M, Pagani F, Yeo K-TJ, et al. Evaluation of imprecision for cardiac troponin assays at low-range concentrations. *Clin Chem* 2004;50(2):327-332.
31. Passing H, Bablok W. A new biometrical procedure for testing the equality of measurements from two different analytical methods. *J Clin Chem Clin Biochem* 1983;21:709-20.
32. US Department of Labor, Occupational Safety and Health Administration, 29 CFR Part 1910.1030, Bloodborne pathogens.
33. US Department of Health and Human Services. *Biosafety in Microbiological and Biomedical Laboratories*. 5th ed. Washington, DC: US Government Printing Office; December 2009.
34. World Health Organization. *Laboratory Biosafety Manual*. 3rd ed. Geneva: World Health Organization; 2004.
35. Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections: Approved Guideline-Fourth Edition*. CLSI Document M29-A4. Wayne, PA: Clinical and Laboratory Standards Institute; 2014.

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©ABBOTT JAPAN LLC 2023