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- 1) Kiesel W. Human plasma Protein C - Isolation, Characterization and Mechanism of Activation by Alpha Thrombin. J. Clin. Invest. 1979; 64:761-769.
- 2) Stenflo J. Structure and Function of Protein C. Seminars in Thrombosis and Haemostasis 1984;10(2): 109-121.
- 3) Francis Jr. RB, Patch MJ. A Functional Assay for Protein C in Human Plasma. Thrombosis Research 1983; 32: 605-613.
- 4) Martinoli JL, Stocker K. Fast Functional Protein C Assay Using Protac, A Novel Protein C Activator. Thrombosis Research 1986; 43: 253-264.
- 5) Clinical and Laboratory Standards Institute. Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Haemostasis Assays; Approved Guideline – Fifth Edition, CLSI Document H21-A5; Vol. 28 No. 5.
- 6) Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.
- 7) Westgard JO, Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process. AACC Press 1986.
- 8) Clinical and Laboratory Standards Institute. Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline. Fourth, Edition, CLSI Document C3-A4; Vol. 26 No.22.
- 9) Griffin JH, Mosher DF, Zimmerman TS, Kleiss AJ. Protein C, an Antithrombotic Protein, is Reduced in Hospitalized Patients with Intravascular Coagulation. Blood 1982; 60: 261-264.
- 10) Mannucci PM, Viganò S. Deficiency of Protein C and Inhibitor of Blood Coagulation. Lancet 1982; 2: 463-467.
- 11) Comp PC, Nixon R, Esmon CT. Determinations of Functional Levels of Protein C, an Antithrombotic Protein, using Thrombin/Thrombomodulin Complex. Blood 1984; 63: 16-21.
- 12) Kottke-Marchant K, Comp P. Laboratory Issues in Diagnosing Abnormalities of Protein C, Thrombomodulin, and Endothelial Cell Protein C Receptor. Arch Pathol Lab med 2002; 126:1337-1348

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