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- 1) Dahlbäck B. Purification of human C4b-binding protein and formation of its complex with vitamin K-dependent protein S, Biochem. J. 1983; 209: 847-856.
- 2) He XH, Shen L, Malmberg AC, Smith KJ. Dahlbäck B, Linse S. Binding site for C4-binding protein in vitamin K-Dependent protein S fully contained in carboxy-terminal laminin-G-type repeats. A study using recombinant factor IX protein S chimeras and surface plasmon resonance. Biochemistry 1997; 36: 3745-3754.
- 3) Serra J, Sales M, Chitolie A, Domènech P, Rossi E, Borrell M, Dahlbäck B. Multicentre Evaluation of IL Test™ Free PS: A Fully Automated Assay to Quantify Free Protein S. Thromb. Haemost. 2002; 88: 975-983.
- 4) Walker FJ. Regulation of Activated Protein C by a New Protein: A Possible Function for Bovine Protein S, J. Biol. Chem. 1980; 255: 5521-5524.
- 5) Suzuki K, Nishioka J. Plasma Protein S Activity measured using Protac, a Snake Venom Derived Activator of Protein C, Thromb. Res. 1988; 49: 214-251.
- 6) Clinical and Laboratory Standards Institute/CLSI. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Hemostasis Assay ; Approved Guideline-Fifth Edition, CLSI/NCCLS Document H21-A5; Vol. 28 No. 5.
- 7) Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.
- 8) Westgard JO, Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press1986.
- 9) Clinical and Laboratory Standards Institute/NCCLS. Preparation and Testing of Reagent Water in the Clinical Laboratory, Third Edition, NCCLS Document C3-A3; Vol. 17 No. 18.
- 10) Liberti G, Bertina RM, Rosendaal FR. Hormonal State rather than Age Influences Cut-off Values of Protein S: Reevaluation of the Thrombotic Risk Associated with Protein S Deficiency. Thromb. Haemost. 1999; 82: 1093-1096.
- 11) Solberg HE. Approved Recommendation (1987) on the Theory of Reference Values. Part 5 Statistical Treatment of Collected Reference Values. Determination of Reference Limits, J. Clin. Chem. Clin. Biochem. 1987; 25: 645-656.
- 12) Sala N, Morell M, Tirado I, Espinosa Y, Llobet D, Fontcuberta J, Soria JM, Volpini V, Estivill X. Linkage disequilibrium between the protein S Heerlen allele and protein S deficiency in Spanish families. Thromb. Haemost. 1995; 73: 1259.
- 13) Maccaferri M, Legnani C, Preda L, Palareti G. Protein S Activity in Patients with Heredofamilial protein S Deficiency and in patients with Juvenile venous Thrombosis. Results of a functional method, Thromb. Res. 1991; 64: 647-658.
- 14) Faioni EM, Valsecchi C, Palla A, Taioli E, Razzari C, Mannucci PM. Free protein S deficiency is a Risk Factor for Venous Thrombosis. Thromb. Haemost. 1997; 78: 1343-1346.
- 15) Mannucci PM, Valsecchi C, Krachmalnicoff A, Faioni EM, Tripodi A. Familial Dysfunction of Protein S, Thromb. Haemost. 1989; 62: 763-766.
- 16) D' Angelo A, Viganò D' Angelo S, Esmon CT, Comp PC. Acquired Deficiencies of Protein S. Protein S Activity during Oral Anticoagulation, in Liver Disease, and in Disseminated Intravascular Coagulation, J.Clin. Invest. 1988; 81: 1445-1454.
- 17) Comp PC, Esmon CT. Recurrent Venous Thromboembolism in Patients with a Partial Deficiency of Protein S, N. Engl. J. Med. 1984; 311: 1522-1528.
- 18) インスツルメンテイション ラボラトリー社 資料
- 19) Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S.Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.

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