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- Gaffney PJ, Creighton LJ, Perry MJ, Callus M, Thorpe R, Spitz M. Monoclonal antibodies to crosslinked fibrin degradation products (XL-FDP) I, Characterization and preliminary evaluation in plasma, British Journal of Haematology. 1988; 68: 83-90.
- Palareti G. Fibrinogen/fibrin Degradation Products: Pathophysiology and Clinical Application, Fibrinolysis. 1993; 7: 60-61.
- Gaffney PJ. The Occurrence and Clinical Relevance of Fibrin Fragments in Blood, Annals New York Academy of Sciences. 1983; 408: 407-423.
- Bounameaux H, de Moerloose P, Perrier A, Reber G. Plasma Measurement of DDimer as Diagnostic Aid in Suspected Venous Thromboembolism: An Overview, Thrombosis and Haemostasis. 1994; 71 (1): 1-6.
- Gaffney PJ, Creighton LJ, Callus M, Thorpe R. Monoclonal antibodies to crosslinked fibrin degradation products (XL-FDP) II, Evaluation in a

variety of clinical conditions. British Journal of Haematology. 1988; 68: 91-96.

- 6) Nolan TE, Smith RP, Devoe LD. Maternal Plasma D-Dimer Levels in Normal and Complicated Pregnancies, Obstetrics & Gynecology. 1993; 81 (2): 235-238.
- 7) Anderson DR, Wells PS. Improvements in the Diagnostic Approach for Patients with Suspected Deep Vein Thrombosis or Pulmonary Embolism, Thrombosis and Haemostasis. 1999; 82 (2): 878-886.
- 8) Kovacs MJ, MacKinnon KM, Anderson D, O' Rourke K, Keeney M, Kearon C, Ginsberg J, Wells PS. A comparison of three rapid D-dimer methods for the diagnosis of Venous Thromboembolism, British J. Haematology. 2001; 115 (1): 140-146.
- 9) Wells PS, Anderson DR, Rodger M, Forgie M, Kearon C, Dreyer J, Kovacs G, Mitchell M, Lewandowski B, Kovacs MJ. Evaluation of D-Dimer in the Diagnosis of Suspected Deep-Vein Thrombosis, N. Engl. J. Med. 2003; 349:1227-35.
- 10) Oswald CT, Menon V, Stouffer GA. The use of D-Dimer in emergency room patients with suspected deep vein thrombosis: a test whose time has come, Journal of Thrombosis and Haemostasis. 2003; 1: 635-636.
- 11) Anderson DR, Kovacs MJ, Kovacs G, Stiell I, Mitchell M, Khoury V, Dryer J, Ward J, Wells PS. Combined use of clinical assessment and D-Dimer to improve the management of patients presenting to the emergency department with suspected deep vein thrombosis (the EDITED study), J Thromb Haemost, 2003; 1: 645-651.
- 12) Curtin N, Highe G, Harris M, Braunstein A, Demattia F, Coss L. Extensive evaluation of the instrumentation laboratory IL test D-Dimer immunoturbidimetric assay on the ACL 9000 determines the D-Dimer cutoff value for reliable exclusion of venous thromboembolism. Lab Hematol. 2004;10 (2): 88-94.
- 13) Legnani C, Pancani C, Palareti G, Guazzaloca G, Coccheri S. Performance of a new, fast D-Dimer test (IL Test D-Dimer) for the management of outpatients with suspected deep vein thrombosis in emergency situations.. Fibrinolysis & Proteolysis 1999. 13(3): 139-141.
- 14) deMoerloose P, Vanrusselt M, Reber G, Arnout J. Performance of the HemosIL Ddimer HS assay for the exclusion of venous thromboembolism, J Thromb Haemost, 2005; 3: 2361-2363.
- 15) Clinical and Laboratory Standards Institute (CLSI). Collection, Transport and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Hemostasis Assays; Approve Guideline, Fifth Edition, CLSI/NCCLS Document H21-A5; Vol. 28 No. 5.
- 16) Kricka LJ. Human Anti-Animal Antibody Interferences in immunological Assays. Clinical Chemistry 1999; 45:7:942-956.
- 17) Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.
- 18) Westgard JO, Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
- 19) Clinical and Laboratory Standards Institute (CLSI). Preparation and Testing of Reagent Water in the Clinical Laboratory; Approved Guideline, Fourth Edition, CLSI/NCCLS Document C3-A4; Vol. 26 No. 22.
- 20) Villa P, Ferrando F, Serra J, Faus H, Mira Y, Vaya A, Aznar J. Quantification of Ddimer using a new fully automated assay: its application for the diagnosis of deep vein thrombosis, Haematologica. 2000 May; 85(5): 520-4.
- 21) Righini M, de Moerloose P, Reber G, Perrier A, Bounameaux H. Should the DDimer Cut-off Value Be Increased in Elderly Patients Suspected of Pulmonary Embolism? Thrombosis and Haemostasis. 2001; 85: 744.
- 22) Hart DJ, Hutchman G, Cuthbert RJG. Evaluation of an automated latex D-Dimmer immunoassay in the clinical assessment of suspected venous thromboembolism.Clin. Lab. Haem. 2002. 24:171-174.. Haem. 2002. 24:171-174.
- 23) Tumanova I, Serra J. New quantitative assay for D-Dimer. Method comparison study: standardization difficulties, Communication presented at the XIII International Congress on Fibrinolysis and Thrombolysis, Barcelona, June 24-28,1996.
- 24) Nieuwenhuizen W. A Reference Material for Harmonisation of D-Dimer Assays, Thrombosis and Haemostasis. 1997; 77 (5): 1031-1033.
- 25) Jennings I, Woods TAL, Kitchen DP, Kitchen S, Walker ID. Laboratory D-Dimmer measurement: Improved agreement between methods through calibration. Thrombosis and Haemostasis. 2007;98; 1127-1135
- 26) Meijer P, Haverkate K, Kluft C, de Moerloose P, Verbruggen B, Spannagl M. A Model for the Harmonisation of Test Results of Different Quantitative D-Dimmer Methods, J Thromb Haemost, 2006; 95: 567-72.
- 27) Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.

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