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- 1. Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.
- 2. Newman DJ, et al. Particle enhanced light scattering immunoassay, Ann Clin Biochem. 1992; 29: 22-42.
- 3. Clinical and Laboratory Standards Institute/NCCLS. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Hemostasis Assays: Approved Guideline - Fifth Edition, CLSI/NCCLS Document H21-A5; Vol. 28 No. 5.
- 4. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.
- 5. Westgard JO, and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.

- 6. Clinical and Laboratory Standards Institute (CLSI). Preparation and Testing of Reagent Water in the Clinical Laboratory: Approved Guideline, Fourth Edition, CLSI document C3-A4; Vol. 26, No. 22.
- 7. Villa P, Fernando F, Serra J, Faus H, Mira Y, Vaya A, Aznar J. Quantification of D-dimer using a new fully automated assay: its application for the diagnosis of deep vein thrombosis. Haematologica. 2000 May; 85(5): 520-4.
- 8. Bounameaux H. Should the D-Dimer Cut-off Value Be increased in Elderly Patients Suspected of Pulmonary Embolism? Thrombosis and Haemostasis. 2001; 85:744.
- 9. Gaffney PJ, et al. Monoclonal antibodies to crosslinked fibrin degradation products (XL-FDP) I, Characterization and preliminary evaluation in plasma, British Journal of Haematology. 1988; 68: 83-90.
- 10. Palareti G. Fibrinogen/fibrin Degradation Products: Pathophysiology and Clinical Application, Fibrinolysis. 1993; 7: 60-61.
- 11. Gaffney PJ. The Occurrence and Clinical Relevance of Fibrin Fragments in Blood, Annals New York Academy of Sciences. 1983; 408: 407-423.
- 12. Bounameaux H, et al. Plasma Measurement of D-Dimer as Diagnostic Aid in Suspected Venous Thromboembolism: An Overview, Thrombosis and Haemostasis. 1994; 71 (1): 1-6.
- 13. Gaffney PJ, et al. Monoclonal antibodies to crosslinked fibrin degradation products (XL-FDP) II, Evaluation in a variety of clinical conditions. British Journal of Haematology. 1988; 68: 91-96.
- 14. Nolan TE, et al. Maternal Plasma D-Dimer Levels in Normal and Complicated Pregnancies, Obstetrics & Gynecology. 1993; 81 (2): 235-238.
- 15. Anderson DR, and Wells PS. Improvements in the Diagnostic Approach for Patients with Suspected Deep Vein Thrombosis or Pulmonary Embolism, Thrombosis and Haemostasis. 1999; 82 (2): 878-886.
- 16. Kovacs MJ, MacKinnon KM, Anderson D, O'Rourke K, Keeney M, Kearon C, Ginsberg J, Wells PS. A comparison of three rapid D-dimer methods for the diagnosis of Venous Thromboembolism, British J. Haematolgy. 2001; 115 (1): 140-146.
- 17. Wells PS, Anderson DR, Rodger M, Forgie M, Kearon C, Dreyer J, Kovacs G, Mitchell M, Lewandowski B and Kovacs MJ. Evaluation of D-Dimer in the Diagnosis of Suspected Deep-Vein Thrombosis, N. Engl. J. Med. 2003; 349:1227-35.
- 18. Oswald CT, Menon V, Stouffer GA. The use of D-Dimer in emergency room patients with suspected deep vein thrombosis: a test whose time has come, Journal of Thrombosis and Haemostasis, 2003; 1: 635-636.
- 19. Anderson DR, Kovacs MJ, Kovacs G, Stiell I, Mitchell M, Khoury V, Dryer J, Ward J and Wells PS. Combined use of clinical assessment and D-Dimer to improve the management of patients presenting to the emergency department with suspected deep vein thrombosis (the EDITED study), J Thromb Haemost. 2003; 1: 645-651.
- 20. Curtin N, Highe G, Harris M, Braunstein A, Demattia F, Coss L. Extensive evaluation of the instrumentation laboratory IL test D-Dimer immunoturbidimetric assay on the ACL 9000 determines the D-Dimer cutoff value for reliable exclusion of venous thromboembolism. Lab Hematol. 2004;10 (2): 88-94.
- 21. Legnani C, Pancani C, Palareti G, Guazzaloca G, Coccheri S. Performance of a new, fast D-Dimer test (IL Test D-Dimer) for the management of outpatients with suspected deep vein thrombosis in emergency situations. Fibrinolysis & Proteolysis 1999; 13(3): 139-141.
- 22. Tumanova I, Serra J. New quantitative assay for D-Dimer. Method comparison study: standardization difficulties, Communication presented at the XIII International Congress on Fibrinolysis and Thrombolysis, Barcelona, June 24-28, 1996.
- 23. Nieuwenhuizen W. A Reference Material for Harmonisation of D-Dimer Assays, Thrombosis and Haemostasis. 1997; 77 (5): 1031-1033.
- 24. Jennings I, Woods TAL, Kitchen DP, Kitchen S, Walker ID. Laboratory D-Dimmer measurement: Improved agreement between methods through calibration. Thromb Haemost. 2007;98: 1127-1135
- 25. Meijer P, Haverkate F, Kluft C, de Moerloose P, Verbruggen B, Spannagl M. A model for the harmonisation of test results of different quantitative D-Dimer methods. Thromb Haemost. 2006; 95: 567-572.

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