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1. Richmond JY, McKinney RW eds. Biosafety in Microbiological and Biomedical Laboratories, U.S. Dept. of Health and Human Services, Public Health Service, 4th Edition, 1999.
2. Newman DJ, et al. Particle enhanced light scattering immunoassay, Ann Clin Biochem. 1992; 29: 22-42.
3. Clinical and Laboratory Standards Institute/NCCLS. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Hemostasis Assays; Approved Guideline - Fifth Edition, CLSI/NCCLS Document H21-A5; Vol. 28 No. 5.
4. Zucker S, Cathey MH, West B. Preparation of Quality Control Specimens for Coagulation, Am. J. Clin. Pathol. 1970; 53: 924-927.
5. Westgard JO, and Barry PL. Cost-Effective Quality Control: Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
6. Clinical and Laboratory Standards Institute. Preparation and Testing of Reagent Water in the Clinical Laboratory, 4th Edition, CLSI Document GP40-A4-AMD.
7. Villa P, Ferrando F, Serra J, Faus H, Mira Y, Vaya A, Aznar J. Quantification of D-dimer using a new fully automated assay: its application for the diagnosis of deep vein thrombosis. Haematologica, 2000 May; 85(5): 520-4.
8. Suspected of Pulmonary Embolism? Thrombosis and Haemostasis. 2001; 85: 744.
9. Gaffney PJ, et al. Monoclonal antibodies to crosslinked fibrin degradation products (XL-FDP) I, Characterization and preliminary evaluation in plasma, British Journal of Haematology. 1988; 68: 83-90.
10. Palareti G. Fibrinogen/fibrin Degradation Products: Pathophysiology and Clinical Application, Fibrinolysis. 1993; 7: 60-61.
11. Gaffney PJ. The Occurrence and Clinical Relevance of Fibrin Fragments in Blood, Annals New York Academy of Sciences. 1983; 408: 407-423.
12. Bounameaux H, et al. Plasma Measurement of D-Dimer as Diagnostic Aid in Suspected Venous Thromboembolism: An Overview, Thrombosis and Haemostasis. 1994; 71 (1): 1-6.
13. Gaffney PJ, et al. Monoclonal antibodies to crosslinked fibrin degradation products (XL-FDP) II, Evaluation in a variety of clinical conditions. British Journal of Haematology. 1988; 68: 91-96.
14. Nolan TE, et al. Maternal Plasma D-Dimer Levels in Normal and Complicated Pregnancies, Obstetrics & Gynecology. 1993; 81 (2): 235-238.
15. Tumanova I, Serra J. New quantitative assay for D-Dimer. Method comparison study: standardization difficulties. Communication presented at the XIII International Congress on Fibrinolysis and Thrombolysis, Barcelona, June 24-28, 1996.
16. Nieuwenhuizen W. A Reference Material for Harmonisation of D-Dimer Assays, Thrombosis and Haemostasis. 1997; 77 (5): 1031-1033.
17. Jennings I, Woods TAL, Kitchen DP, Kitchen S, Walker ID. Laboratory D-Dimer measurement: Improved agreement between methods through calibration. Thromb Haemost. 2007;98: 1127-1135.
18. Meijer P, Haverkate F, Kluft C, de Moerloose P, Verbruggen B, Spannagl M. A model for the harmonisation of test results of different quantitative D-Dimer methods. Thromb Haemost, 2006; 95: 567-572.

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