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- 1) Henry JB. Clinical diagnosis and management by laboratory methods. Philadelphia: WB Saunders,1979.
- 2) Clinical and Laboratory Standards Institute. Collection, Transport, and Processing of Blood Specimens for Testing Plasma-Based Coagulation and Molecular Hemostasis Assays; Approved Guideline - Fifth Edition, CLSI Document H21-A5; Vol. 28 No. 5
- 3) Clinical and laboratory Standards Institute. One Stage Prothrombin Time (PT) Test and Activated Partial Thromboplastin Time (APTT) Test; Approved Guideline - Second Edition, CLSI Document H47-A2; Vol. 28 No. 20.
- 4) Zucker S, Cathey MH, West B. Preparation of Quality Control Specimen for Coagulation, Am J Clin Pathol 1970; 53: 924-927.
- 5) Westgard JO, Barry PL. Cost-Effective Quality Control; Managing the Quality and Productivity of Analytical Process, AACC Press 1986.
- 6) Sirridge MS. Laboratory evaluation of hemostasis. Philadelphia: Lea & Febiger,1967.
- 7) Shapiro GA, Huntzinger SW, Wilson JE. Variation among commercial activated partial thromboplastin time reagents in response to heparin. Am J Clin Pathol 1977; 67: 477.
- 8) Ratnoff OD, Crum JD. Activation of Hageman factor by solutions of ellagic acid. J Lab Clin Med 1964; 63: 359.
- 9) Van den Besselaar AMHP, Neuteboom J, Bertina RM. Effect of Synthetic Phospholipids on the response of the Activated Partial Thromboplastin Time to Heparin. Blood Coag Fibrin 1993; 4: 895-903.
- 10) Hirsh J, Raschke R, Warkentin TE, Dalen JE, Deykin D, Poller L. Heparin: mechanism of action, pharmacokinetics, dosing considerations, monitoring, efficacy and safety. Chest 1995 Oct; 108 (4 Suppl): 258S-275S.
- 11) Brill-Edwards P, Ginsburg JS, Johnson M, Hirsch J. Establishing a Therapeutic Range for Heparin Therapy. Ann Int Med 1993; 119: 104-109.
- 12) CLSI (Clinical and Laboratory Standards Institute) EP05-A3. Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline-Third Edition

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