LDHI

093319/03-2601

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1. Burtis CA, Bruns DE *Tietz Fundamentals of Clinical Chemistry and Molecular Diagnostics*. ed. 7. Philadelphia: WB Saunders; 331-332; 2015.
2. Quist J, Hill AR. Serum lactate dehydrogenase (LDH) in Pneumocystis carinii pneumonia, tuberculosis, and bacterial pneumonia. *Chest*. 1995; 108:415-418.
3. Winzelberg GG, Hull JD, Agar JW, Rose BD, Pletka PG. Elevation of serum lactate dehydrogenase levels in renal infarction. *JAMA* 1979; 242: 268-269.
4. Zein JG et al., Prognostic Significance of Elevated Serum Lactate Dehydrogenase (LDH) in Patients with Severe Sepsis; *Chest*. 2004; 126(4_MeetingAbstracts) :873S
5. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37 °C, Part 3. Reference procedure for the measurement of catalytic concentration of lactate dehydrogenase. *Clinical Chemistry and Laboratory Medicine*, 40(6), 643-648. 2002.
6. Doumas BT, et al. Differences Between Values for Plasma and Serum in Tests Performed in the Ektachem 700 XR Analyzer, and Evaluation of "Plasma Separator Tubes (PST)." *Clin. Chem.* 35: 151-153; 1989.
7. Peake MJ, et al. Mechanism of Platelet Interference with Measurement of Lactate Dehydrogenase Activity in Plasma. *Clin. Chem.* 30:518-520; 1984.
8. Greenberg N, Byrne D. Plasma Lactate Dehydrogenase Activity Assayed with the Kodak Ektachem 700 Analyzer Is Unaffected by Platelet Contamination. *Clin. Chem.* 31: 1022; 1985.
9. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests. Ed. 2. Washington D.C.: AACC Press; 3-335; 1997.
10. Clinical Laboratory Handbook for Patient Preparation and Specimen Handling. Fascicle VI: Chemistry/Clinical Microscopy. Northfield, IL: College of American Pathologists; 1992.
- * 11. Tietz NW. *Tietz Textbook of Clinical Chemistry*. ed. 3. Philadelphia: WB Saunders; 671; 1999.
- * 12. Calam RR. Specimen Processing Separator Gels: An Update. *J Clin Immunoassay*. 11:86-90; 1988.
- * 13. Volz M, Kellner H-M. Kinetics and Metabolism of Pyrazolones (Propyphenazone, Aminopyrine and Dipyrone). *Br. J. Clin Pharmacol.* 10: 299S-308S; 1980.
- * 14. Eldor A, Zylber-Katz E, Levy M. The Effect of Oral Administration of Dipyrone on the Capacity of Blood Platelets to Synthesize Thromboxane A2 in Man. *Eur. J. Clin. Pharmacol.* 26: 171-176; 1984.
- * 15. Bagnoud M-A, Reymond J-Ph. Interference of Metamizol (Dipyrone) on the Determination of Creatinine with the Kodak Dry Chemistry Slide Comparison with the Enzymatic Method from Boehringer. *Eur. J. Clin. Chem. Clin. Biochem.* 31: 753-757; 1993.
- * 16. Comandone A, Passera R High dose methotrexate in adult patients with osteosarcoma: Clinical and pharmacokinetic results, *Acta Oncologica*, 44:4, 406-411; 20053.
17. Young DS. Effects of Drugs on Clinical Laboratory Tests. ed. 4. Washington D. C.: AACC Press; 1995.
18. Friedman RB, Young DS. Effects of Disease on Clinical Laboratory Tests. Washington, D.C.: AACC Press; 1990.
- * 19. Siemens Healthcare Diagnostics, Inc. Lactate Dehydrogenase L-P (LDLP) for ADVIA® Chemistry Systems Instructions for Use. Siemens Healthcare Diagnostics Inc., Tarrytown, NY 10591-5097 USA. Publication 04346317 Rev. C, 2008-10.

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