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1. Aminotransferases in Diseases, Robert Rej, Clinics in Laboratory Medicine. Vol 9, No 4, December 1989.
2. Evaluation of abnormal liver-enzyme results in asymptomatic patients, Daniel S. Pratt and Marshall M. Kaplan, N Engl J Med 2000; 342:1266-1271, April 27, 2000.
3. Tietz NW (ed). Fundamentals of Clinical Chemistry. ed. 3. Philadelphia: WB Saunders; 369-371; 1987.
4. Schumann et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Part 4. Reference procedure for the measurement of catalytic concentrations of alanine aminotransferase. Clin. Chem. Lab. Med. 40(7): 718-724; 2002.
5. Schumann et al. IFCC primary reference procedures for the measurement of catalytic activity concentrations of enzymes at 37°C. International Federation of Clinical Chemistry and Laboratory Medicine (IFCC). Part 5. Reference procedure for the measurement of catalytic concentrations of aspartate aminotransferase. Clin. Chem. Lab. Med. 40(7):725-733; 2002.
6. Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests. ed. 5. Washington D.C.: AACC Press; 3-94;2000.
7. Calam RR. Specimen Processing Separator Gels: An Update. J Clin Immunoassay. 11:86-90; 1988.
8. Young DS. Effects of Drugs on Clinical Laboratory Tests. ed. 4. Washington D.C.: AACC Press; 1995.
9. Friedman RB, Young DS. Effects of Disease on Clinical Laboratory Tests. Washington, D.C.: AACC Press; 1990.
10. Lippi, G. et al. Influence of hemolysis on routine clinical chemistry testing. Clinical Chemical Laboratory Medicine, 44(3), pp. 311-316; 2006.
11. Thomas, L., Muller, M., Schumann, G., Klein, G., Lunau, S., Pick, K.-H., & Oswald, S. (2005). Consensus of DGKL and VDGH for interim reference intervals on enzymes in serum. Journal of Laboratory Medicine, 29(5), 301-308.

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