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1. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. Sixth edition. Edited by N Rafai, AR Horvath, CT Wittwer. Elsevier, 2018.
2. Pearson TA, Mensah GA, Alexander RW, et al. Markers of inflammation and cardiovascular disease: application to clinical and public health practice: A statement for healthcare professionals from the Centers for Disease Control and Prevention and the American Heart Association. Circulation. 2003; 107 (3) :499-511.
3. Bard RL, Rubenfire M, Eagle K, Clarke NS, Brook RD. Utility of C-reactive protein measurement in risk stratification during primary cardiovascular disease prevention. Am J Cardiol. Jun 1; 95 (11) :1378-1379; 2005.
4. Ridker PM, Stampfer MJ, Rifai N. Novel Risk Factors for Systemic Atherosclerosis: A Comparison of C-Reactive Protein, Fibrinogen, Homocysteine, Lipoprotein (a), and Standard Cholesterol Screening as Predictors of Peripheral Arterial Disease. JAMA. 2001; 285 (19) :2481-2485.
5. Heesch C, Hamm CV, Bruehmer J, Simoons ML. Predictive value of C-reactive protein and troponin T in patients with unstable angina: a comparative analysis. J Am Coll Cardiol. 2000 May; 35 (6) :1535-42.
6. Ridker PM, Rifai N, Pfeffer MA, et al. Inflammation, pravastatin, and the risk of coronary events after myocardial infarction in patients with average cholesterol levels. Cholesterol and Recurrent Events (CARE) Investigators. Circulation. 1998; 98:839-844.
7. Bowen, RAR and Remaley, AT. Interferences from blood collection tube components on clinical chemistry assays. Biochem Med. 24:31-44; 2014.
8. Guder WG, et al. Quality of Diagnostic Samples. Recommendations of the Working Group on Preanalytical Quality of the German Society for Clinical Chemistry and Laboratory Medicine. Darmstadt, Germany; 2010; pages 42-43.
9. 2. WHO/DIL/LAB/99.1 Rev.2 World Health Organization. Use of Anticoagulants in Diagnostic Laboratory Investigations 2002; page 28.
10. Young, D.S. Effects of Drugs on Clinical Laboratory Tests. Ed. 5. Washington D.C.: AACC Press; 2000
11. Friedman, R.B.; Young, D.S. Effects of Disease on Clinical Laboratory Tests. Ed. 4. Washington, D.C.: AACC Press; 2001.
12. Levinson SS. The Nature of Heterophilic Antibodies and Their Role in Immunoassay Interference. J Clin Immunoassay. 15:108-115; 1992.
13. Dati F, Johnson AM, Whicher JT. The existing interim consensus reference ranges and the future approach. Clin Chem Lab Med. 2001; 39 (11) :1134-6.
14. H. Schlebusch, N. Liappis, E. Kalina, C. Klein: High Sensitive CRP and Creatinine: Reference Intervals from Infancy to Childhood. Laboratoriums Medizin. 2002; 26 (5) p.341-346.
15. European Commission, The Certification of a Matrix Reference Material for Immunochemical Measurement of 15 Serum Proteins, ERM-DA470, Report EUR 15243 EN and 16882 EN, European Communities; 2004.
16. CLSI. Interference Testing in Clinical Chemistry. 3rd ed. CLSI Guideline EP07. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.

17. CLSI. Supplemental Tables for Interference Testing in Clinical Chemistry. 1st ed. CLSI supplement EP37. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.

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