A1C1

095207/02-2601

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1. A. Finke, U. Kobold, W. Hoelzel, C. Weykamp, K. Miedema, J-O. Jeppsson. Approved IFCC Reference Method for the Measurement of HbA1c in Human Blood. Clin Chem Lab Med; 40(1): 78-89, 2002.
2. Goldstein DE, Little RR, Wiedmeyer HM, England JD, Rohlfing CL, Wilke AL. Is glycohemoglobin testing useful in diabetes mellitus: Lessons from the diabetes control and complications trial. Clin Chem; 40:1637-1640, 1994.
3. The Diabetes Control and Complications Trial Research Group: The effect of intensive treatment of diabetes on the development and progression of long-term complications in insulin-dependent diabetes mellitus. N Engl J Med 329:977-86; 1993.
4. Intensive blood glucose control with sulphonylureas or insulin compared with conventional treatment and risk of complications in patients with type 2 diabetes (UKPDS 33) . UK Prospective Diabetes Study (UKPDS) Group. Lancet 325:837-53; 1998.
5. 糖尿病治療ガイド2022-2023、日本糖尿病学会
6. Bowen, RAR and Remaley, AT. Interferences from blood collection tube components on clinical chemistry assays. Biochem Med. 24:31-44; 2014.
7. CLSI. Supplemental Tables for Interference Testing in Clinical Chemistry. 1st ed. CLSI supplement EP37. Wayne, PA: Clinical and Laboratory Standards Institute; 2018.
8. Kashiwagi A., et al. International clinical harmonization of glycated hemoglobin in Japan: From Japan Diabetes Society to national glycohemoglobin standardization program values. Diabetology International. 3:8-10.
9. Kashiwagi A., et al. International clinical harmonization of glycated hemoglobin in Japan: From Japan Diabetes Society to national glycohemoglobin standardization program values. Journal of Diabetes Investigation. 3:39-40.
10. <https://www.niddk.nih.gov/health-information/professionals/clinical-tools-patient-management/diabetes/sickle-cell-trait-hemoglobinopathies-diabetes>. Accessed on August 20, 2021.
11. <http://www.ngsp.org/factors.asp>. Accessed on August 20, 2021.
12. International Association of Diabetes and Pregnancy Study Groups Recommendations on the Diagnosis and Classification of Hyperglycemia in Pregnancy. Diabetes Care. 33:No 3; 2010.
13. International Expert Committee Report on the Role of the A1c Assay in the Diagnosis of Diabetes. Diabetes Care. 32: (7): 1327-1334; 2009.
14. Standards of Medical Care in Diabetes. Diabetes Care. 38: Supplement 1; S8-S16, 2015.
15. http://www.who.int/diabetes/publications/report-hba1c_2011.pdf. Accessed on August 20, 2021.

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