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- 1) Landsteiner K. Über agglutinationserscheinungen normalen menschlichen blutes. Klin Wschr 14: 1132: 1901.
- 2) Moersch C. Neue tatsachen uber die blutkorperchen agglutinationen. Zbl Bakt 46: 49, 456: 1908.
- 3) Coombs RRA, Mourant AE, Race RR. A new test for the detection of weak and incomplete Rh agglutinins. Br J Exp Pathol 26: 255: 1945.
- 4) Levine P. Cross-matching technics, especially the Coombs test. Proc NY State Assn of Public Health Laboratories 35: 9: 1955.
- 5) Dacie JV, Crookston JH, Christenson WN, Incomplete cold antibodiesrole of complement in sensitization to antiglobulin serum by potentially haemolytic antibodies. Br J Haematol 3: 77: 1957.
- 6) Hill G, Pollack W, Levine P. A critical evaluation of a rapid 'broad spectrum compatibility test. Proc Seventh Congress Int Soc Blood Trans, Rome, 1958.
- 7) Pollack W, Hager HJ, Hollenbeck LL Jr, The specificity of anti-human gamma globulin reagents. Transfusion 2: 17: 1962.
- 8) Stroup M, MacIlroy M.Evaluation of the albumin antiglobulin technic in antibody detection. Transfusion 5: 184: 1965.
- 9) Blood group antigens & antibodies as applied to hemolytic disease of the newborn.Raritan, NJ: Ortho Diagnostic Systems, 1968.
- 10) Austen KF et al. Nomenclature of complement. Immunochemistry 7: 137: 1970.
- 11) Issitt PD, Smith TR, Evaluation of antiglobulin reagents. In: A seminar on performance evaluation.Washington, DC: American Association of Blood Banks, 61, 65-6: 1976.
- 12) Garratty G. Laboratory investigation of drug-induced immune hemolytic anemia and/or positive direct antiglobulin tests.Washington, DC: American Association of Blood Banks, 1980.
- 13) Howard JE, Winn LC, Grumet FC, Garratty G, Petz LD. Clinical significance of the anti-complement component of antiglobulin antisera. Transfusion 22: 269: 1982.
- 14) Technical manual. 9th ed. Arlington, VA: American Association of Blood Banks, 91-93: 1985.
- 15) Kosmin M, Treacy M. Clinical significance of the direct antiglobulin test.Raritan, NJ: Ortho Diagnostic Systems, 1978.

- 16) Giles C, Engelfriet CP. Working party on the standardization of antiglobulin reagents of the expert panel of serology. Vox Sang 38: 178-9: 1980.
- 17) Moore JA, Chaplin H Jr. Anti-C3d antiglobulin reagents. II: Preparation of an antiglobulin serum monospecific for C3d, Transfusion 14: 416: 1974.
- 18) Köhler G, Milstein C, Continuous cultures of fused cells secreting antibody of predefined specificity, Nature 256: 495: 1975.
- 19) SANDRA J. Nance, GEORGE Garratty A new potentiator of Red Blood Cell Antigen-Antibody reactions, A. J. C. P. May 1987.

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