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- 1.平成24年 人口動態統計（厚生労働省大臣官房統計情報部）
- 2.Wang MC, Valenzuela LA, Murphy GP and Chu TM. Purification of a human prostate specific antigen. Invest Urol. 17:159-163; 1979.
- 3.Wang MC, Papsidero LD, Kuriyama M et al. Prostate antigen: a new potential marker for prostatic cancer. The Prostate.2:89-96; 1981.
- 4.Papsidero LD, Wang MC, Valenzuela LA et al. A Prostate Antigen in Sera of Prostatic Cancer Patients. Cancer Res.40:2428-2432; 1980.
- 5.Oesterling JE. Prostate specific antigen: a critical assessment of the most useful tumor marker for adenocarcinoma of the prostate. J Urol. 145:907-923; 1991.
- 6.Stamey TA, Kabalin JN, McNeal JE et al. Prostate specific antigen in the diagnosis and treatment of adenocarcinoma of the prostate. II. Radical prostatectomy treated patients. J Urol. 141:1076-1083; 1989.
- 7.Chybowski FM, Larson Keller JJ, Bergstralh EJ and Oestriling JE. Predicting radionuclide bone scan findings in patients with newly diagnosed, untreated prostate cancer: prostate specific antigen is superior to all other clinical parameters. J Urol. 145:313-318; 1991.
- 8.Stamey TA, Kabalin JN, Ferrari M. Prostate specific antigen in the diagnosis and treatment of adenocarcinoma of the prostate. III. Radiation treated patients. J Urol. 141:1084-1087; 1989.

- 9.Stamey TA, Kabalin JN, Ferrari M, Yang N. Prostate specific antigen in the diagnosis and treatment of adenocarcinoma of the prostate. IV. Anti-androgen treated patients. J Urol. 141:1088-1090; 1989.
- 10.Catalona WJ, Partin AW, Slawin KM, Brawer MK, Flanigan RC, Patel A, et al.. Use of the percentage of free prostate specific antigen to enhance differentiation of prostate cancer from benign prostatic disease: A prospective multicenter clinical trial. JAMA. 279: 1542-1547;1998.
- 11.Potter, SR, Reckwitz, T, and Partin AW. The use of percent free PSA for early detection of prostate cancer. J. Andrology. 20/4: 449-453;1999.
- 12.Catalona WJ, Smith DS, Wolfert RL, Wang TJ, Rittenhouse HG, Ratliff TL, et al. Evaluation of percentage of free serum prostate specific antigen to improve specificity of prostate cancer screening. JAMA. 274:1214-1220; 1995.
- 13.Prestigiacoma AF, Liua H, Pettersson, K, Wolfert RL, Stamey TA. A Comparison of the Free Fraction of Serum Prostate Specific Antigen in men with Benign and Cancerous Prostates: The Best Case Scenario. The Journal of Urology 1996; 156:350-354.
- 14.Summers M et al. Luminogenic Reagent Using 3-Chloro 4-Hydroxy Acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. Clin Chem. 41:S73; 1995.
- 15.Kumari, G.R., and Malati, T. (2004) Stability of Total and Free Prostate Specific Antigen in Serum Samples at Different Storage Conditions. Indian Journal of Clinical Biochemistry. 19 (2), 10-13.
- 16.Sokoll, L.J., Bruzek, D.J., Dua, R., Dunn, W., Mohr, P., Wallerson, G., Eisenberger, M, Partin, A.W. and Chan, D.W. (2002) Short-Term Stability of the Molecular Forms of Prostate-Specific Antigen and Effect on Percent Complexed Prostate-Specific Antigen and Percent Free Prostate-Specific Antigen. Urology. 60 (Supplement 4A), 24-30.
- 17.Calam RR. Specimen Processing Separator Gels: An Update. J Clin Immunoassay. 11:86-90; 1988
- 18.Levinson, SS.: The Nature of Heterophilic Antibodies and Their Role in Immunoassay Interference. J Clin Immunoassay 15: 108-115 (1992).
- \* 19.Clinical and Laboratory Standards Institute (CLSI). Interference Testing in Clinical Chemistry; Approved Guideline, Second Edition. CLSI document EP07-A2. Wayne, PA: CLSI; 2005.
- \* 20.Livanou E, et.al.: Serum biotin levels in patients undergoing chronic hemodialysis. Nephron. 1987; 46(3):331-2.
- \* 21.Paul Grimsey et.al.: Population pharmacokinetics of exogenous biotin and the relationship between biotin serum levels and in vitro immunoassay interference. International Journal of Pharmacokinetics. 2017; 2 (4): 247-56.
- \* 22.European Commission Scientific Committee on Food: Opinion of the Scientific Committee on Food on the Tolerable Upper Intake Level of Biotin (expressed on 26 September 2001). Rue de la Loi 200, B-1049 Bruxelles/Wetstraat 200, B-1049 Brussels-Belgium.
- \* 23.Mock and Mock: Serum concentrations of bisnorbiotin and biotin sulfoxide increase during both acute and chronic biotin supplementation J Lab Clin Med 1997;129:384-8.
- \* 24.Marcelo C. Batista: Biotin interference in immunoassays mimicking subclinical Graves' disease and hyperestrogenism: a case series. Clin Chem Lab Med. 2017 May 1;55(6):e99-e103.
- 25.Thiel RP, Oesterling JE, Wojno KJ, Partin AW, Chan DW, Carter HB, et al. A multicenter comparison of the diagnostic performance of free prostate-specific antigen. Urology 1996;48(6A):45-50.
- 26.Young DS. Effects of Drugs on Clinical Laboratory Tests. ed. 4. Washington, D.C.: AACC Press; 1995.
- 27.Friedman RB, Young DS. Effects of Disease on Clinical Laboratory Tests. ed. 3. Washington, D.C.: AACC Press; 1997.
- 28.Tryding N, Tufvesson C, Sonntag O (eds). Drug Effects in Clinical Chemistry. ed. 7. Stockholm: The National Corporation of Swedish Pharmacies, Pharmasoft AB, Swedish Society for Clinical Chemistry; 1996.
- 29.Guess, H. A. et al. : J. Urol., 155(1) : 3, 1996
- 30.Andriole GL et al : J Urol 175(5) : 1657-1662, 2006
- 31.Marks LS et al : J Urol 176(3) : 868-874, 2006
- 32.Passing HW, A New Biometrical Procedure of Testing the Equality of Measurements from Two Different Analytical Methods. J Clin Chem Biochem. 21:709-720 1983.

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