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- Summers M. et al.: Luminogenic Reagent Using 3-Chloro 4-Hydroxy Acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. Clinical Chemistry. 41. S73; 1995.
- Calam RR. Specimen Processing Separator Gels: An Update. J Clin Immunoassay. 11:86-90; 1988.
- Levinson SS. The Nature of Heterophilic Antibodies and Their Role in Immunoassay Interference. J Clin Immunoassay. 15:108-115; 1992.
- Müller B, Becker KL, Schächinger H, et al. Calcitonin precursors are reliable markers of sepsis in a medical intensive care unit. Crit Care Med. 28 (4) :977-983; 2000.
- Harbarth S, Holeckova K, Froidevaux C, et al. Diagnostic value of procalcitonin, interleukin-6 and interleukin-8 in critically ill patients admitted with suspected sepsis. Am J Respir Crit Care Med. 164:396-402; 2001.
- Brunkhorst FM, Wegscheider K, Forycki ZF, et al. Procalcitonin for early diagnosis and differentiation of SIRS, sepsis, severe sepsis, and septic shock. Intensive Care Med. 26:148-152; 2000.
- American College of Chest Physicians / Society of Critical Care Medicine Consensus Conference Committee. Definitions for sepsis and organ failure and guidelines for the use of innovative therapies in sepsis. Crit Care Med. 20 (6) :864-874; 1992.
- Meisner M. Procalcitonin - Biochemistry and Clinical Diagnosis. UNIMED Science. ISBN 978-3-8374-1241-3; 2010.
- Reith H.B, et al. Procalcitonin in early detection of postoperative complications. Dig Surg. 15 (3) :260-5; 1998.
- Meisner M, et al. Postoperative plasma concentrations of procalcitonin after different types of surgery. Intensive Care Med. 24 (7) :680-4; 1998.
- Chiesa C. et al. Reliability of procalcitonin concentrations for the diagnosis of sepsis in critically ill neonates. Clinical Infectious Diseases. 26:664-72; 1998.
- Kruger S. and Welte T. Biomarkers in community-acquired pneumonia. Expert Rev Respir Med. 6 (2) :203-14; 2012.
- Vijayan AL, Maya V, Ravindran S, Saikant R, Lakshmi S, Kartik R and Manoj G. Procalcitonin: a promising diagnostic marker for sepsis and antibiotic therapy. Journal of Intensive Care. (5) :51: 1 - 7; 2017.
- Markus B, Peter AW. The inflammatory response in sepsis. Trends Immunol. 34 (3) :129-36; 2013.
- Lever A, Mackenzie I. Sepsis: definition, epidemiology, and diagnosis. Br Med J. 335:879-83; 2007.
- Dunja M, Snezana B, Arsen U, Biljana D, Vladimir V. Use of presepsin and procalcitonin for prediction of SeptiFast results in critically ill patients. J Crit Care. 40:197-201; 2017.
- Angus DC, Van der Poll T. Severe sepsis and septic shock. N Engl J Med. 369:840-51; 2013.
- Moore LJ, McKinley BA, Turner K. The epidemiology of sepsis in general surgery patients. J Trauma. 70 (3) :672-80; 2011.
- Taylor R, Jones A, Kelly S, et al. A review of the value of procalcitonin as a marker of infection. Cureus. Apr 10;9 (4) :e1148. doi:10.7759/cureus.1148; 2017.
- Becker KL, Nylén ES, White JC, Müller B, Snider RH. Clinical review 167: Procalcitonin and the calcitonin gene family of peptides in inflammation, infection, and sepsis: a journey from calcitonin back to its precursors. J Clin Endocrinol Metab. 89: 1512-1525; 2004.
- Karzai W, Oberhoffer M, Meier-Hellmann A, Reinhart K. Procalcitonin - a new indicator of the systemic response to severe infections. Infection. 25: 329-334; 1997.
- Shieraw B, Ebisa B, Kumar K, Boutin A, Frieri M. The Role of Procalcitonin as a Biomarker in Sepsis. J Infect Dis Epidemiol. 2:006; 2016.
- Morgenthaler N.G, Struck J. et al. Detection of procalcitonin (PCT) in healthy controls and patients with local infection by a sensitive ILMA. Clin Lab. 48 (5-6) :263-270; 2002.
- Christ-Crain M, Jaccard-Stolz D. et al. Effect of procalcitonin-guided treatment on antibiotic use and outcome in lower respiratory tract infections: cluster-randomised, single-blinded intervention trial. Lancet. 363 (9409) :600-607; 2004.
- Luyt C.E. et al. Procalcitonin kinetics as a prognostic marker of ventilator-associated pneumonia. Am. J. Resp. Crit. Care Med. 171 (1) 48-53; 2005.
- Brunkhorst F.M. et al. Kinetics of procalcitonin in iatrogenic sepsis. Intensive Care Med. 24: 888-892; 1998.
- Schuetz P, Christ-Crain M. et al. Effect of procalcitonin-based guidelines vs standard guidelines on antibiotic use in lower respiratory tract infections: the ProHOSP randomized controlled trial. JAMA. 302 (10) :1059-1066; 2009.
- Bouadma L, Luyt C.E. et al. Use of procalcitonin to reduce patients' exposure to antibiotics in intensive care units (PRORATA trial) : a multicentre randomised controlled trial. Lancet. 375 (9713) :463-474; 2010.
- Schuetz P. et al. Serial Procalcitonin Predicts Mortality in Severe Sepsis Patients: Results From the Multicenter Procalcitonin Monitoring SEpsis (MOSES) Study. Crit. Care Med. 45 (5) :781-789; 2017.
- Deming WE, Statistical Adjustment of Data. New York, NY: John Wiley and Sons; 1943.

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