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- Summers M et al. Luminogenic Reagent Using 3-Chloro 4-Hydroxy Acetanilide to Enhance Peroxidase/Luminol Chemiluminescence. Clin Chem. 41:S73; 1995.
- Calam RR. Specimen Processing Separator Gels: An Update. J Clin Immunoassay. 11:86-90; 1988.
- Young DS. Effects of Drugs on Clinical Laboratory Tests. ed. 4. Washington, D.C.: AACC Press, 1995.
- Friedman RB, Young DS. Effects of Disease on Clinical Laboratory Tests. ed. 3. Washington, D.C.: AACC Press; 1997.
- Tryding N, Tufvesson C, Sonntag O (eds). Drug Effects in Clinical Chemistry. Ed. 7. Stockholm: The National Corporation of Swedish Pharmacies, Pharmasoft AB, Swedish Society for Clinical Chemistry; 1996.
- Robert-Guroff M, Weiss SH, Giron JA, et al. Prevalence of antibodies to HTLV-I, -II, and -III in intravenous drug abusers from an AIDS endemic region. JAMA 1986;255(2):3133-3137.
- Lee HH, Swanson P, Rosenblatt JD, et al. Relative prevalence and risk factors of HTLV-I and HTLV-II infection in US blood donors. Lancet 1991;337:1435-1439.
- Okochi K, Sato H, Hinuma Y. A retrospective study on transmission of adult T-cell leukemia virus by blood transfusion: seroconversion in recipients. Vox Sang 1984;46:245-253.

- Sato H, Okochi K. Transmission of human T-cell leukemia virus (HTLV-I) by blood transfusion: demonstration of proviral DNA in recipients' blood lymphocytes. Int J Cancer 1986;37:395-400.
- Kaplan JE, Litchfield B, Rouault C, et al. HTLV-I associated myelopathy associated with blood transfusion in the United States: epidemiologic and molecular evidence linking donor and recipient. Neurology 1991;41:192-197.
- Saxton EH, Lee H, Swanson P, et al. Detection of human T-cell leukemia/lymphoma virus type I in a transfusion recipient with chronic myelopathy. Neurology 1989;39:841-844.
- Donegan E, Pell P, Lee H, Shaw GM, Mosley JW. Transmission of human T-lymphotropic virus type I by blood components from a donor lacking anti-p24: a case report. The Transfusion Safety Study Group. Transfusion 1992;32:68-71.
- Hino S, Yamaguchi K, Katimine S, et al. Mother-to-child transmission of human T-cell leukemia virus type I. Jpn J Cancer Res 1985;76:474-480.
- Kajiya W, Kashiwagi S, Ikematsu H, Hayashi J, Nomura H, Okochi K. Intrafamilial transmission of adult T cell leukemia virus. J Infect Dis 1986;154:851-857.
- Sullivan MT, Williams AE, Fang CT, Notari EP, Poiesz BJ, Ehrlich GD. Human T-lymphotropic virus (HTLV) types I and II infection in sexual contacts and family members of blood donors who are seropositive for HTLV type I or II. Transfusion 1993;33:585-590.
- Centers for Disease Control (CDC). Human T-lymphotropic virus type I screening in volunteer blood donors — United States, 1989. MMWR Morb Mortal Wkly Rep 1990;39:915-924.
- Blattner WA, Takatsuki K, Gallo RC. Human T-cell leukemia-lymphoma virus and adult T-cell leukemia. JAMA 1983;250:1074-1080.
- Hinuma Y, Komoda H, Chosa T, et al. Antibodies to adult T-cell leukemia-virus-associated antigen (ATLA) in sera from patients with ATL and controls in Japan: a nationwide sero-epidemiologic study. Int J Cancer 1982;29:631-635.
- Gessain A, Barin F, Vermant JC, et al. Antibodies to human T-lymphotropic virus type-I in patients with tropical spastic paraparesis. Lancet 1985;2:407-410.
- Rodgers-Johnson P, Gajdusek DC, Morgan OSC, Zaninovic V, Sarin PS, Graham DS. HTLV-I and HTLV-II antibodies and tropical spastic paraparesis. Lancet 1985;2:1247-1248.
- Osame M, Usuku K, Izumo S, et al. HTLV-I associated myelopathy, a new clinical entity. Lancet 1986;1:1031-1032.
- Gessain A, Gout O. Chronic myelopathy associated with human T-lymphotropic virus type I (HTLV-I). Ann Intern Med 1992;117:933-946.
- Wiley CA, Nerenberg M, Cros D, Soto-Aguilar MC. HTLV-I polymyositis in a patient also infected with the human immunodeficiency virus. N Engl J Med 1989;320:992-995.
- Blattner WA. Epidemiology of HTLV-I and associated diseases. In: Blattner WA, ed. Human retrovirology: HTLV. New York: Raven Press, 1990:251-66.
- Newton M, Cruickshank K, Miller D, et al. Antibody to human T-lymphotropic virus type-I in West-Indian-born UK residents with spastic paraparesis. Lancet 1987;1:415-416.
- The T-and B-cell malignancy study group. Statistical analyses of clinico-pathological, virological and epidemiological data on lymphoid malignancies with special references to adult T-cell leukemia/lymphoma: a report of the second nationwide study of Japan. Jpn J Clin Oncol 1985;15:517-535.
- Saxinger W, Blattner WA, Levine PH, et al. Human T-cell leukemia virus (HTLV-I) antibodies in Africa. Science 1984;225:1473-1476.
- Maloney EM, Biggar RJ, Neel JV, Taylor ME, et al. Endemic human T-cell lymphotropic virus type II infection among isolated Brazilian Amerindians. J Infect Dis 1992;166:100-107.
- Babona DV, Nurse GT. HTLV-I antibodies in Papua New Guinea. Lancet 1988;2:1148-1149.
- Yanagihara R, Jenkins CL, Alexander SS, Mora CA, Garruto RM. Human T lymphotropic virus type I infection in Papua New Guinea: high prevalence among the Hagahai confirmed by western analysis. J Infect Dis 1990;162:649-654.
- Bastian I, Gardner J, Webb D, Gardner I. Isolation of a human T-lymphotropic virus type I strain from Australian aboriginals. J Virol 1993;67:843-851.
- Kirkland MA, Frasca J, Bastian I. Adult T-cell leukaemia; lymphoma in an aborigine. Aust NZ J Med 1991;21:739-741.
- Rosenblatt JD, Golde DW, Wachsman W, et al. A second isolate of HTLV-II associated with atypical hairy-cell leukemia. New Engl J Med 1986;315:372-377.

34. Murphy EL, Engstrom JW, Miller K, Sacher RA, Busch MP, Hollingsworth CG. HTLV-II associated myelopathy in a 43-year-old woman. REDS Investigators. Lancet 1993;341:757-758.
35. Jacobson S, Lehky T, Nishimura M, Robinson S, McFarlin DE, Dhib-Jalbut S. Isolation of HTLV-II from a patient with chronic progressive neurological disease clinically indistinguishable from HTLV-I-associated myelopathy/tropical spastic paraparesis. Ann Neurol 1993;33:392-396.
36. Hjelle B, Appenzeller O, Mills R, Alexander S et al. Chronic neurodegenerative disease associated with HTLV-II infection. Lancet 1992;339:645-646.
37. Hjelle B, Mills R, Swenson S, Mertz G, Key C, Allen S. Incidence of hairy cell leukemia, mycosis fungoides, and chronic lymphocytic leukemia in first known HTLV-II-endemic population. J Infect Dis 1991;163:435-440.
38. Maloney EM, Biggar RJ, Neel JV, et al. Endemic human T-cell lymphotropic virus II infection among isolated Brazilian Amerindians. J Infect Dis 1992;166:100-107.
39. Wiktor SZ, Cannon RD, Atkinson WL, et al. Infection with human T-lymphotropic virus types I and II in sexually transmitted disease clinics in Baltimore and New Orleans. J Infect Dis 1992;165:920-924.
40. Zella D, Mori L, Sala M, et al. HTLV-II infection in Italian drug abusers. Lancet 1990;336:575-576.
41. Hjelle B, Wilson C, Cyrus S, et al. Human T-cell leukemia virus type II infection frequently goes undetected in contemporary US blood donors. Blood 1993;81:1641-1644.

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