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1. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition. NCCLS document H3-A5 [ISBN 1-56238-515-1]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
2. Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. NCCLS, 940 West Valley Road, Suite1400, Wayne, PA 19087-1898 USA, 2004.
3. Tietz, NW, Clinical Guide to Laboratory Tests -- 3rd ed. Philadelphia: W. B. Saunders Co., 1995: 186-188.
4. Clinical and Laboratory Standards Institute/NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard -- Fifth Edition; NCCLS document H1-A5 [ISBN 1-56238-519-4]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
5. Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline – Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
6. National Cholesterol Education Program. Third report of the National Cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation and Treatment of High Blood Cholesterol in Adults (ATP III) Final Report. NIH Publication No. 02-5215, 2002.
7. Stadtman, TC. Methods in Enzymology, Vol. III, Colowick, SP, and Caplan, NO, (Eds.), Academy Press, New York, NY, 1957: pp 392–394, 678–681.
8. Flegg, HM. An investigation of the determination of serum cholesterol by an enzymatic method. Ann Clin Biochem 1973; 10:79–84.
9. Roschlau P, Bernt EM, Gruber W. Enzymatische Bestimmung des Gesamt-Cholesterins im Serum. Z Klin. Chem. Klin. Biochem. 1974; 12: 403–407.
10. Rautela, GS, Liedtke, RJ. Automated enzymic measurement of total cholesterol in serum. Clin Chem 1978; 24.
11. Tietz NW, Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Burtis CA, Ashwood ER, Bruns DE, Elsevier Saunders, 2006; pp 903-943,1695.
12. Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline- Second Edition. CLSI/NCCLS Document EP5-A2 (ISBN 1- 56238-542-9) Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA, 19087, 2004.
13. Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline- Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.

**14. Wu, Alan H.B. Tietz Clinical Guide to Laboratory Tests, Fourth Edition, 1995, WB Saunders Company, St. Louis, MO, 1995, p. 244.

**15. Rifai, N, Warnick GR, Dominiczak MH. Handbook of Lipoprotein Testing, Second Edition, 2000, American Association for Clinical Chemistry, Washington DC, p. 176.

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