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- Tietz, NW. Clinical Guide to Laboratory Tests, 3rd ed. Philadelphia; W.B. Saunders Company, 1995: p. 374.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Collection　of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Fifth Edition.　NCCLS document H3-A5 [ISBN 1-56238-515-1]. NCCLS,940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 2nd ed. Washington, DC: AACCPress, 1997: p. 4-499.
- Clinical and Laboratory Standards Institute /NCCLS. Tubes and Additives for Venous Blood Specimen Collection; Approved Standard—Fifth Edition. NCCLS document H1-A5 [ISBN 1-56238-519-4]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Kroll, M, Evaluating Interference Caused by Lipemia, Clin Chem 2004; 50: 1968-1969.
- Kaplan LA, Psece AJ. Clinical Chemistry Theory, Analysis, and Correlation, 3rd ed. St. Louis;Mosby, Inc., 1996: p 699, 713-714.
- Burtis, CA, Ashwood ER, Bruns DE. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed. St. Louis: Elsevier Saunders, 2005, 1186-1193.
- Burtis CA, Ashwood ER. Tietz Fundamentals of Clinical Chemistry, 5th ed. Philadelphia: W.B. Saunders Company, 2001: pp. 31-41 (specimen collection), 597-598.
- Smith, FE, Herbert, J, Gaudin, J, Hennessy, J, Reid, GR, Serum iron determination using ferene triazine, Clin Biochem 1984; 17:306–310.
- Higgins, T, Novel chromogen for serum iron determinations, Clin Chem 1981; 27:1619–1620.
- Artiss, JD, Vinogradov, S, Zak, B, Spectrophotometric study of several sensitive reagents for serum iron, Clin Biochem 1981; 14:311–315.
- Artiss, JD, Standbergh, DR, Zak, B, Study of continuous flow automation for serum iron on comparing several sensitive reagents, Microchem J 1983; 28:275–284.
- Hennessy, DG, Reid, GR, Smith FE, Thompson, SL, Ferene—a new spectrophotometric reagent for iron, Can J Chem 1984; 62:721–724.
- Clinical and Laboratory Standards Institute/NCCLS. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline- Second Edition. CLSI/NCCLS Document EP5-A2 (ISBN 1- 56238-542-9) Clinical and Laboratory Standards Institute, 940 West Valley Road, Suite 1400, Wayne, PA, 19087, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Second Edition. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. NCCLS, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002

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