SIEMENS

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- Burtis, CA, Ashwood ER. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed, Philadelphia; W.B. Saunders Co., 2006; 41-45.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard-Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute /NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard-Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Friedman RB, Young DS. Effects of Disease on Clinical Laboratory Tests, 3rd ed, Washington DC; AACC Press, 1997; 3-228-229.
- Berth M, Delanghe, J. Protein Precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of the literature, Acta clinica Belgica 2004; 59-5.
- Smogorzewska A, Flood JG, Long WH, Dighe AS. Paraprotein interference in automated chemistry analyzers. Clin Chem 2004;50(9) 1691-3.
- Clinical and Laboratory Standards Institute/NCCLS. *Interference Testing in Clinical Chemistry; Approved Guideline-Second Edition*. CLSI/NCCLS document EP7-A2 [ISBN 1-56238-584-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2005.
- Rifai N, Warnick GR, Dominiczak MH. Handbook of lipoprotein testing, Washington: AACC Press 1997; 131 (methodology), 13 (cholesterol transport).
- Genest J. Genetics and prevention: a new look at high density lipoprotein cholesterol. Cardiol. In Rev. 2002;10(1):61-71.
- Genest J, Frohlich J, Fodor G, McPherson R. Recommendations for the management of dyslipidemia and the prevention of cardiovascular disease: 2003 update, CMAJ 2003; 169 (9).
- Fletcher B, et al. AHA Scientific Statement Managing Abnormal Blood Lipids: A collaborative approach, Circulation 2005; 112:3184 - 3209.
- Third report of the National Cholesterol Education Program (NCEP) expert panel on detection, evaluation, and treatment of high blood cholesterol in adults (Adult Treatment Panel III) final report. Circ. 2002;106; 3232, 3172.
- Clinical and Laboratory Standards Institute /NCCLS. *Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition*. CLSI/NCCLS document EP5-A2 [ISBN 1-56238-542-9]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Clinical and Laboratory Standards Institute/NCCLS. *Method Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Second Edition*. CLSI/NCCLS document EP9-A2 [ISBN 1-56238-472-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2002.
- Clinical and Laboratory Standards Institute/NCCLS. *Protocols for Determination of Limits of Detection and Limits of Quantitation; Approved Guideline*. CLSI/NCCLS document EP17-A [ISBN 1-56238-551-8]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- http://www.cdc.gov/labstandards/crmln.htm.

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