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- Yamanishi H. et al. Fully automated measurement of total iron-binding capacity in serum. Clin Chem 43:12, 2413-2417;1997.
- Higgins T. Novel chromogen for serum iron determinations. Clin Chem 1981;27:1619-1620.
- Artiss JD, Vinogradov S, Zak B. Spectrophotometric study of several sensitive reagents for serum iron, Clin Biochem 1981;14:311-315.
- Hennessy DG, Reid GR, Smith FE, Thompson SL. Ferene—a new spectrophotometric reagent for iron. Can J Chem 1984;62:721-724.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard – Fifth Edition*. CLSI/NCCLS document H3-A5 [ISBN 1-56238-515-1]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Tubes and Additives for Venous Blood Specimen Collection; Approved Standard – Fifth Edition*. CLSI/NCCLS document H1-A5 [ISBN 1-56238-519-4]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2003.
- Clinical and Laboratory Standards Institute/NCCLS. *Procedures for the Handling and Processing of Blood Specimens; Approved Guideline – Third Edition*. CLSI/NCCLS document H18-A3 [ISBN 1-56238-555-0]. CLSI, 940 West Valley Road, Suite 1400, Wayne, PA 19087-1898 USA, 2004.
- Burtis CA, Ashwood ER. Tietz Textbook of Clinical Chemistry, W.B. Saunders Co., Philadelphia, PA 1999, Third Edition pp 54-56.
- Young DS. Effects of Preanalytical Variables on Clinical Laboratory Tests, 3rd ed. Washington, DC: AACC Press, 2007;p 570.
- Clinical Practice Guidelines for the Treatment of Anemia of Chronic Renal Failure; Guideline 5. National Kidney Foundation, Inc., 1997.
- Tietz NW. Clinical Guide to Laboratory Tests W.B.Saunders Co., Philadelphia, PA 1990.
- Tietz NW. Textbook of Clinical Chemistry and Molecular Diagnostics, 4th ed., Burtis CA, Ashwood ER, Bruns DE, Elsevier Saunders, 2006; pp 1186-1193.

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