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- Fiske CH, Subbarow U. *The colorimetric determination of phosphorus*, J Biol Chem 1925; 66:375.
- Daly, JA, Ertlingshausen, G, *Direct Method for Determining Inorganic Phosphate in Serum with the “Centrifichem”*, Clin Chem 1972;18:263-265.
- CLSI. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition*. CLSI document GP41-A6 [ISBN 1-56238-650-6] Wayne, PA: Clinical and Laboratory Standards Institute; 2007.
- CLSI. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition*. CLSI document GP39-A6 [ISBN 1-56238-740-5]. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- CLSI. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition*. CLSI document GP44-A4 [ISBN 1-56238-724-3]. Wayne, PA: Clinical and Laboratory Standards Institute; 2010.
- Kaplan LA, Pesce AJ. *Clinical Chemistry Theory Analysis and Correlation*, Third edition, Mosby Co, St. Louis, MO, 1996; p 552.
- CLSI. *Urinalysis; Approved Guidelines—Third Edition*. CLSI document GP16-A3 [ISBN 1-56238-687-5]. Wayne, PA: Clinical and Laboratory Standards Institute; 2009.
- CLSI. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition*. CLSI document EP07-A2 [ISBN: 1-56238-584-4] Wayne, PA: Clinical and Laboratory Standards Institute; 2005.
- Tietz NW. Ed. *Clinical Guide to Laboratory Tests*, 4th Ed. WB Saunders Co., Philadelphia PA, 2006:852-853.
- CLSI. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guidelines—Third Edition*. CLSI document EP28-A3c [ISBN: 1-56238-682-4] Wayne, PA: Clinical and Laboratory Standards Institute, 2010.
- Most S, Pollard G, Van Wormer D, Steige H, Lembke L, and Jones JD. An Evaluation of a Phosphorus Method for the Du Pont aca®, Du Pont Company, Wilmington, DE, September, 1976.
- Burtis CA, Ashwood ER, Bruns DE. *Tietz Text book of Clinical Chemistry and Molecular Diagnostics*, 5th edition, Elsevier Inc. Philadelphia, PA, 2012: pp 145 - 162 (Specimen Collection and Processing), p 2163 (Reference Intervals and Values).
- CLSI. *Defining, Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guidelines—Second Edition*. CLSI document EP05-A2 [ISBN: 1-56238-542-9] Wayne, PA: Clinical and Laboratory Standards Institute, 2004.
- CLSI. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guidelines—Third Edition*. CLSI document EP09-A3 [ISBN: 1-56238-888-6] Wayne, PA: Clinical and Laboratory Standards Institute, 2013.
- CLSI. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guidelines—Second Edition*. CLSI document EP17-A2 [ISBN: 1-56238-796-0] Wayne, PA: Clinical and Laboratory Standards Institute, 2013.
- Burtis CA, Ashwood ER, Bruns DE. Tietz Text book of Clinical Chemistry and Molecular Diagnostics, 5th edition, Elsevier Inc. Philadelphia, PA, 2012: pp 145 - 162 (Specimen Collection and Processing), p 1751 (Measurement of Phosphorus), p 2163 (Reference Intervals and Values).
- Tauskey, H.H., Shorr, E., J.Biol Chem 1953, 202:675-685.
- Guder, WG, Narayanan S, Wisser H, et al. List of analytes – Preanalytical variables. Annex In: Samples: From the patient to the Laboratory. Darmstadt, Germany: GIT Verlag; 1996: Annex 20 – 1, 38 – 9.

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