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1. Pruden E.L., Siggaard-Andersen O., and Tietz N.W., Chapter 30 (Blood Gases and pH), of 『Tietz Textbook of Clinical Chemistry, Second Edition』, ed. C.A.Burtis and E.R.Ashwood.W.B.Saunders Company, Philadelphia, 1994.
2. Gambino et al., 「The measurement of CO₂ content with the autoanalyzer」, Am.J.Clin Path, 45, p. 406,

- 1966.
3. Ungerer et al., [Discordance between measured and calculated total carbon dioxide] , Clin.Chem., 36 (12), p. 2093-2096, 1990.
4. D.B.Sacks, Chapter 25 (p. 837) of [Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, fourth edition] , C.A.Burtis, E.R.Ashwood, and D.E.Bruns eds., Elsevier Saunders, St. Louis, 2006.
5. Chapter 141, Blood Glucose of J. Michael McMillin, Walker HK, Hall WD, Hurst JW, editors.Clinical Methods: The History, Physical, and Laboratory Examinations, third edition Boston: Butterworths; 1990.
6. M.G. Scott, V.A. LeGrys, J.S. Klutts, Chapter 27 (p.985) of Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, fourth edition, C.A. Burtis, E.R. Ashwood, D.E. Bruns eds, Elsevier Saunders, St. Louis, 2006.
7. D. Young, E. Bermes, Jr, Chapter 2 of Tietz Textbook of Clinical Chemistry, third edition, C.A. Burtis and E.R. Ashwood, eds., W.B. Saunders Company, Philadelphia, 1999.
8. CLSI. [Procedures for the Collection of Arterial Blood Specimens; Approved Standard] , H11-A4 (ISBN 1-56238-545-3), 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2004.
9. CLSI. [Interference Testing in Clinical Chemistry; Approved Guideline] , EP07-A2 (ISBN 1-56238-480-5), 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2002.
10. M.G.Scott et al, [Electrolytes and Blood Gases] , Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, 4th Ed, Ed.C.A.Burtis, E.R.Ashwood, D.E.Bruns, Page 991, 2006.
11. D.B.Endres and R.K.Rude, Chapter 49 (p. 1901) of [Tietz Textbook of Clinical Chemistry and Molecular Diagnostics, fourth edition] , C.A.Burtis, E.R.Ashwood, and D.E.Bruns eds., Elsevier Saunders, St. Louis, 2006.
12. Whillier S., Raftos J.E., Chapman B., Kuchel P.W. (2009). [Role of N-acetylcysteine and cystine in glutathione synthesis in human erythrocytes] .Redox Report: Communications in Free Radical Research, 14 (3), 115.
13. Ventura P., Panini R., Pasini M.C., Scarpetta G., Salvioli G. (1999). [N-Acetyl Cysteine Reduces Homocysteine Plasma Levels After Single Intravenous Administration by Increasing Thiols Urinary Excretion] . Pharmacological Research.40 (4), 345-350.
14. CLSI. [Interference Testing in Clinical Chemistry; Approved Guideline] , EP07-A2 (ISBN 1-56238-480-5), 940 West Valley Road, Suite 1400, Wayne, Pennsylvania 19087-1898 USA, 2002.
15. Richie J.P., Nichenametla S., Neidig W., Calcagnotto A., Haley J.S., Schell T.D., Muscat J.E. (2015).Randomized controlled trial of oral glutathione supplementation on body stores of glutathione, Eur J Nutr., 54: 251. B9. Allen J. and Bradley R.D. (2011).Effects of Oral Glutathione Supplementation on Systemic Oxidative Stress Biomarkers in Human Volunteers, J Altern Complement Med.; 17(9): 827-833.
16. Allen J. and Bradley R.D. (2011).Effects of Oral Glutathione Supplementation on Systemic Oxidative Stress Biomarkers in Human Volunteers, J Altern Complement Med.; 17(9): 827-833.
17. Guerci B, B. M. (2003).Accuracy of an electrochemical sensor for measuring capillary blood ketones by fingerstick samples during metabolic deterioration after continuous subcutaneous insulin infusion interruption in type 1 diabetic patients.Diabetes Care., 26(4), 1137-1141.
18. Laffel L. (1999).Ketone bodies: A review of physiology, pathophysiology and application of monitoring to diabetes., Diabetes Metab Res Rev., 15(6), 412-426.
19. Chiu RW., Ho CS., Tong SF., Ng KF., Lam CW. (2002) Evaluation of a new handheld biosensor for point-of-care testing of whole blood beta-hydroxybutyrate concentration., Hong Kong Med J., 8(3), 172-176.
20. LaGow B et al., eds.PDR Lab Advisor.A Comprehensive Point-of-Care Guide for Over 600 Lab Tests. First ed. Montvale, NJ: Thomson PDR; 2007.
21. Bristol-Myers Squibb. Hydrea (hydroxyurea capsules, USP) prescribing information. Princeton, NJ; March 2001.
22. Rodriguez, G.I.; Kuhn, J.G.; Weiss, G.R.; Hilsenbeck, S.G.; Eckardt, J.R.; Thurman, A.; Rinaldi, D.A.; Hodges, S.; Von Hoff, D.D.; Rowinsky, E.K. (1998).A bioavailability and pharmacokinetic study of oral and intravenous hydroxyurea.Blood, 91(5), 1533-1541.
23. Whillier S., Raftos J.E., Chapman B., Kuchel P.W. (2009). [Role of N-acetylcysteine and cystine in glutathione synthesis in human erythrocytes] .Redox Report: Communications in Free Radical Research, 14(3), 115.
24. Ventura P., Panini R., Pasini M.C., Scarpetta G., Salvioli G. (1999). [N-Acetyl Cysteine Reduces Homocysteine Plasma Levels After Single Intravenous Administration by Increasing Thiols Urinary Excretion] . Pharmacological Research.40(4), 345-350.
25. S.M.Wendrot, T.N.Heady, D.M.Haverstick, L.M.Bachmann, M.G.Scott, J.C. Boyd, and D.E.Bruns. Falsely increased chloride and missed anion gap elevation during treatment with sodium thiosulfate. Clinica Chimica Acta 2014; 431: 77-79.
26. B.E.Statland, [Clinical Decision Levels for Lab Tests, Medical Economic Books] , Oradell, NJ, 1987.
27. Pruden E.L., Siggaard-Andersen O., and Tietz N.W., Chapter 30 (Blood Gases and pH), of [Tietz Textbook of Clinical Chemistry, Second Edition] , ed. C.A.Burtis and E.R.Ashwood.W.B.Saunders Company, Philadelphia, 1994.

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