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1. Clinical and Laboratory Standards Institute. Procedures for the Handling and Processing of Blood Specimens; Approved Guideline—Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. CLSI Document H18-A3.
2. Data on file at Siemens Healthcare Diagnostics.
3. Clinical and Laboratory Standards Institute. Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP7-A2.
4. Holick MF. Vitamin D Deficiency. N Engl J Med. 2007;357:266–81.
5. Holick MF. MrOs is D-ficient. J Clin Endocrinol Metab. 2009;94(4):1092-3.
6. Rollins G. Vitamin D Testing—What’s the Right Answer? Labs Grapple with Confusing Analytics, Evidence. Clinical Laboratory News. July 2009;35(7): 1,6.
7. Freeman R. Vitamin D: The sunshine hormone. How and when to treat deficiencies. Menopausal Medicine. May 2009; S8-11.
8. Misra M, Pacaud D, Petrak Z, Collett-Solberg PF, Kappy M. Vitamin D Deficiency in Children and Its Management: Review of Current Knowledge and Recommendations, Pediatrics. 2008; 122:398-417
9. Holick MF. The role of vitamin D for bone health and fracture prevention. Curr Osteoporosis Rep. 2006;4(3):96-102.

10. Okazaki R, Ozono K, Fukumoto S, Inoue D, Yamauchi M, Minagawa M, Michigami T, Takeuchi Y, Matsumoto T, Sugimoto T. Assessment criteria for vitamin D deficiency/in sufficiency in Japan – proposal by an expert panel supported by Research Program of Intractable Diseases, Ministry of Health, Labour and Welfare, Japan, The Japanese Society for Bone and Mineral Research and The Japan Endocrine Society. Endocrine Journal Advance Publication.
11. Özkan B, Hatun, Bereket A. Vitamin D intoxication. Turk J Pediatr. 2012;54:93-98.
12. Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
13. Boscatto LM, Stuart MC. Heterophilic antibodies: a problem for all immunoassays. Clin Chem. 1988;34:27-33.
14. Inloes R, Clark, D and Drobnies A: Interference of fluorescein, used in retinal angiography with certain clinical laboratory tests. Clin Chem 1987, 33:2126-2127.
15. Blutt SE, Weigel NL. Vitamin D and prostate cancer. Proc Soc Exp Biol Med. 1999; 221(2):89-98.
16. Feskanich D, Ma J, Fuchs CS, et al. Plasma vitamin D metabolites and risk of colorectal cancer in women. Cancer Epidemiol Biomarkers Prev. 2004; 13(9):1502-1508.
17. Shin MH, Holmes MD, Hankinson SE, Wu K, Colditz GA, Willett WC. Intake of dairy products, calcium, and vitamin D and risk of breast cancer. J Natl Cancer Inst. 2002; 94(17):1301-1311.
18. Krause R, Buhning M, Hopfenmuller W, Holick MF, Sharma AM. Ultraviolet B and blood pressure. Lancet. 1998; 352(9129):709-710.
19. Lips P. Vitamin D deficiency and secondary hyperparathyroidism in the elderly: consequences for bone loss and fractures and therapeutic implications. Endocr Rev. 2001; 22(4):477-501.
20. Lips P, Hosking D, Lippuner K, et al. The prevalence of vitamin D inadequacy amongst women with osteoporosis: an international epidemiological investigation. J Unrestricted Med. 2006; 260(3):245-254.
21. Hypponen E, Laara E, Reunanen A, Jarvelin MR, Virtanen SM. Intake of vitamin D and risk of type 1 diabetes: a birth-cohort study. Lancet. 2001; 358(9292): 1500-1503.
22. Munger KL, Zhang SM, O' Reilly E, et al. Vitamin D intake and incidence of multiple sclerosis. Neurology. 2004; 62(1):60-65.
23. Merlino LA, Curtis J, Mikuls TR, Cerhan JR, Criswell LA, Saag KG. Vitamin D intake is inversely associated with rheumatoid arthritis: results from the Iowa Women' s Health Study. Arthritis Rheum. 2004; 50(1):72-77.
24. Hollis BW, Horst RL. The Assessment of Circulating 25(OH)D and 1,25(OH)2D: Where We Are and Where We Are Going. J Steroid Biochem Mol Biol. March 2007;103(3-5):473-76.
25. Hart GR. Medical Conditions Associated with Vitamin D Deficiency and the Clinical Consequences. Review Series Volume 1, Sept 2004;1-8.
26. National Institutes of Health - Office of Dietary Supplements. OSD Vitamin D Initiative -Vitamin D Standardization Program (VDSP). Accessed February 2, 2015. <http://ods.od.nih.gov/Research/VitaminD.aspx#vdsp>
27. Sempos CT, Vesper HW, Phinney KW, Thienpont LM, Coates PM, The Vitamin D Standardization Program(VDSP). Vitamin D Status as an International Issue: National Surveys and the Problem of Standardization. Scandinavian Journal of Clinical & Laboratory Investigation, 2012;72(Suppl 243):32-40.
28. Thienpont LM, Stepan HCM, Vesper HW. Standardization of Measurements of 25-Hydroxyvitamin D3 and D2. Scandinavian Journal of Clinical & Laboratory Investigation, 2012; 72(Suppl 243): 41-49.
29. Centers for Disease Control and Prevention Standardization of Measurement Procedures. Accessed February 2, 2015. http://www.cdc.gov/labstandards/hs_standardization.html
30. Clinical and Laboratory Standards Institute. Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline. Wayne, PA: Clinical and Laboratory Standards Institute; 2003. CLSI Document EP6-A.
31. Clinical and Laboratory Standards Institute. Evaluation of Precision Performance of Quantitative Measurement Methods; Approved Guideline-Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2004. CLSI Document EP5-A2.
32. Clinical and Laboratory Standards Institute. Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline-Second Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.
33. Centers for Disease Control. 1988. Update: Universal precautions for prevention of transmission of human immunodeficiency virus, hepatitis B virus and other bloodborne pathogens in healthcare settings. MMWR, 37:377-382, 387, 388.
34. Clinical and Laboratory Standards Institute. Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline-Third Edition. Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document M29-A3.
35. Federal Occupational Safety and Health Administration, Bloodborne Pathogens Standard, 29 CFR 1910.1030.
36. Clinical and Laboratory Standard Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline-Third Edition*. Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.

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