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- Panteghini M, et al. Diagnostic value of four assays for lipase determination in serum: A comparative reevaluation. *Clin Biochem.* 1991;24:497-503.
- Young DS. *Effects of Preanalytical Variables on Clinical Laboratory Tests.* 2nd ed. Washington, DC: AACC Press; 1997:3–354.
- Clinical and Laboratory Standards Institute. *Protection of Laboratory Workers From Occupationally Acquired Infections; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document M29-A4.
- Clinical and Laboratory Standards Institute. *Procedures for the Collection of Diagnostic Blood Specimens by Venipuncture; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2007. CLSI Document GP41-A6.
- Clinical and Laboratory Standards Institute. *Tubes and Additives for Venous and Capillary Blood Specimen Collection; Approved Standard—Sixth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP39-A6.
- Burtis CA, Ashwood ER. *Tietz Textbook of Clinical Chemistry and Molecular Diagnostics.* 4th ed. Philadelphia, PA: W.B. Saunders Co.; 2006:41–56.
- Clinical and Laboratory Standards Institute. *Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests; Approved Guideline—Fourth Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document GP44-A4.
- Guder W, Fonseca-Wolheim W, Heil O et al. Maximum permissible transport and storage times for analysis of blood (serum, plasma) urine and cerebrospinal fluid. *DG Klinische Chemische Mitteilungen* 1995;26:207-224.
- Zhang DJ, Elswick RK, Miller WG, Bailey JL. Effect of serum-clot contact time on clinical chemistry results. *Clin Chem.* 1998 (44):1325-1333.
- Wu AHB, ed. *Tietz Clinical Guide to Laboratory Tests.* 4th ed. St. Louis, MO: WB Saunders Company;2006:676.
- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP07-ED3.

- Clinical and Laboratory Standards Institute. *Interference Testing in Clinical Chemistry; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2005. CLSI Document EP07-A2.
- Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Data on file at Siemens Healthcare Diagnostics.
- Young DS. *Effects of Drugs on Clinical Laboratory Tests.* 3rd ed. Washington, DC: AACC Press; 1990.
- Smogorzewska A, Flood JG, Long WH, Dighe AS. *Paraprotein interference in automated chemistry analyzers.* *Clin Chem* 2004;50(9) 1691-3.
- Berth M, Delanghe, J. *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins:* 2 case reports and a review of the literature, *Acta clinica Belgica* 2004;59(5):263-273.
- Whitcomb DC, Lowe ME. Human pancreatic digestive enzymes. *Dig Dis Sci.* 2007 Jan;52(1):1-17.
- Lott JA, Patel ST, Sawhney AK, Kazmierczak SC, Love, JE. Assays of serum lipase: analytical and clinical considerations. *Clin Chem.* 1986; 32/7:1290-1302.
- Neumann U, et. al. A sensitive colorimetric assay for the kinetic lipase determination in serum (Boehringer Mannheim Chemicals). Abstract 13th Int. *Congress for Clin. Chem.* (ICCC) 28.6.-3.7. 1987 Den Haag, Netherlands.
- Clinical and Laboratory Standards Institute. *Evaluation of Precision of Quantitative Measurement Procedures; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2014. CLSI Document EP05-A3.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2013. CLSI Document EP09-A3.
- Clinical and Laboratory Standards Institute. *Defining, Establishing, and Verifying Reference Intervals in the Clinical Laboratory; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2010. CLSI Document EP28-A3c.
- Clinical and Laboratory Standards Institute. *Evaluation of the Linearity of Quantitative Measurement Procedures: A Statistical Approach; Approved Guideline.* Wayne, PA: Clinical and Laboratory Standards Institute; 2003. CLSI Document EP06-A.
- Clinical and Laboratory Standards Institute. *Measurement Procedure Comparison and Bias Estimation Using Patient Samples; Approved Guideline—Third Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2018. CLSI Document EP09c.
- Clinical and Laboratory Standards Institute. *Evaluation of Detection Capability for Clinical Laboratory Measurement Procedures; Approved Guideline—Second Edition.* Wayne, PA: Clinical and Laboratory Standards Institute; 2012. CLSI Document EP17-A2.

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